

# ELBOW AND FOREARM ORTHOSES

Available in adult and pediatric sizes.



Forearm Pronation



Forearm Supination



Glove acts as a web space splint.

## COMFORT COOL™ SPIRAL ARM SPLINT



*Provides thumb abduction and dynamic forearm pronation or supination.*

- The 2" (5.1cm) wide spiral forearm strap is made of 1/8" (3.2mm) loop neoprene and easily attaches to the glove using Velcro® hook.
- Helps decrease abnormal tone and increase joint range of motion in persons with neurological disorders, orthopedic conditions and brachial plexus injuries.
- The strap supplies the dynamic force needed for forearm supination or pronation, depending on its spiral placement on the arm. Neoprene can be trimmed for a custom fit. Applies a prolonged, gentle pull that places the upper extremity into a more functional position.
- Includes low-cut glove with wraparound strap, which assists with thumb abduction. This **Web Space Glove** is also available separately as a replacement, or it can be used on its own.
- Ideal for widening thumb web space due to muscle tone imbalance or contractures. Glove is made of 1/16" (1.6mm) perforated, loop neoprene with a terry cloth liner for comfort and coolness.
- Hand wash and air dry.
- Latex free.

Left	Right	MP Circumference
NC34800	NC34801	Pediatric 5 1/2" to 6 1/4" (14 to 16cm)
NC34802	NC34803	Small 6 1/4" to 7" (16 to 18cm)
NC34804	NC34805	Medium 7" to 8" (18 to 20cm)
NC34806	NC34807	Large 8" to 9" (20 to 23cm)

### Web Space Glove Replacement

Left	Right	
NC34800-H	NC34801-H	Pediatric
NC34802-H	NC34803-H	Small
NC34804-H	NC34805-H	Medium
NC34806-H	NC34807-H	Large



Removable thermoplastic insert limits flexion and allows elbow extension.

## COMFORT COOL™ ULNAR NERVE ELBOW PROTECTOR WITH GEL PAD

*Allows greater functional use of the arm because the elbow is not totally immobilized.*

- Made of 1/16" (1.6mm) perforated neoprene with a soft terrycloth liner. Circumferential design with open proximal end is easy to put on and adjust.
- Provides protection and support for bursitis, ulnar nerve irritation, tendinitis and post-surgical repair.
- Mold the precut, 1/8" (3.2mm) Prism™ splint insert to custom fit the elbow at the flexion limit. Design allows the elbow to move from full extension to the limited flexion position, protecting the ulnar nerve from damaging tension.
- Removable gel pad protects the vulnerable posterior elbow from painful bumps and pressure points.
- Hand wash and air dry.
- Latex free.
- If between sizes, choose the larger size.



*Includes thermoplastic insert, neoprene support and gel elbow pad.*

		Forearm Circumference
NC33850-1	Small	8" to 9 1/2" (20 to 24cm)
NC33850-2	Medium	9 1/2" to 11" (24 to 28cm)
NC33850-3	Large	11" to 12 1/2" (28 to 32cm)
NC33850-4	X-Large	Over 12 1/2" (Over 32cm)



## PIL-O-SPLINT® ELBOW SUPPORT

*Soft, adjustable wrap promotes comfortable elbow extension.*

- Made of breathable foam covered with soft fabric for cushioned support.
- Helps relieve pain associated with cubital tunnel syndrome, sprains, tendinitis and tennis elbow.
- Dual Velcro® straps allow easy compression adjustments with one hand.
- The removable, hard plastic stay lies lateral to the elbow and helps block elbow flexion.
- **Standard** measures 11 1/4" (29cm) long.
- **X-Large** measures 13" (33cm) long.
- Machine washable.
- Latex free.



Removable plastic stay.

		Forearm Circumference
NC10113	Standard	Up to 11" (28cm)
NC10115	X-Large	Over 11" (28cm)



## HEELBO® ELBOW PROTECTOR

*Conforms around bony prominences.*

- **Standard** model has an air-cushioned polyurethane foam pad with a tricot lining.
- **Gel Pad** model provides additional protection against shock and vibration.
- Protects elbow or heel from abrasion and promotes healing of pressure ulcers.
- Machine or hand wash. Air dry.
- Sold in pairs. Both styles are latex free.

Standard	Gel Pad		Elbow/Heel Circumference
NC15961-1	NC15961	Small	7" (18cm)
NC15961-2	NC15962	Medium	8" (20cm)
NC15961-3	NC15963	Large	8½" (22cm)
NC15961-4	NC15964	X-Large	9½" (24cm)
NC15961-6	NC15966	XX-Large	11½" (29cm)
Standard Heelbo™ with Foam Pad		(Pair)	
Heelbo™ with Gel Pad		(Pair)	

*Now available in larger sizes.*



## NORCO™ ELBOW/HEEL PROTECTOR

*Economical elbow and heel protection.*

- Available in two padding styles that protect the elbow and/or heel from pressure sores or skin ulceration due to shearing forces.
- Protector offers excellent skin aeration and is ideal for intermittent use.
- The **Foam Elbow/Heel Protector** has a sewn-in foam pad that measures 6" (15cm) in diameter.
- The **Gel Elbow/Heel Protector** has a removable, shock-absorbing gel pad that measures 3½" x 6¼" (8.9 x 16cm).
- Sleeves are machine washable. Remove gel pad prior to washing. Air dry.
- Sold in pairs. Both styles are latex free.
- To size, measure circumference at the elbow crease or heel. If in between sizes, choose the larger size.

Foam Pad		Circumference Stretches to
NC58801-1	Small	Less 9" (Less 23cm)
NC58801-2	Med	9" to 11" (23cm to 28cm)
NC58801-3	Large	11" to 13" (28cm to 33cm)
NC58801-4	X-Large	13" to 15" (33cm to 38cm)
NC58801-5	XX-Large	15" to 18" (38cm to 46cm)

Gel Pad		Circumference Stretches to
NC58800-1	Small	Less 9" (Less 23cm)
NC58800-2	Med	9" to 11" (23cm to 28cm)
NC58800-3	Large	11" to 13" (28cm to 33cm)
NC58800-4	X-Large	13" to 15" (33cm to 38cm)
NC58800-5	XX-Large	15" to 18" (38cm to 46cm)



## NEOPRENE ELBOW SUPPORT

*Protects the elbow from hard surfaces.*

- Neoprene support measures 8" (20cm) long. Pad measures ¼" (13mm) thick.
- Use for treating pain or pressure caused by bursitis or nerve inflammation.
- Comfortable sleeve provides warmth to help reduce pain. Contoured neoprene pad provides support and shock absorption.
- Hand or machine wash. Air dry.
- Latex free.
- To size, measure circumference at the widest part of the forearm.

		Forearm Circumference
NC15433-0	X-Small	8" to 9" (20 to 23cm)
NC15433-1	Small	9" to 10" (23 to 25cm)
NC15433-2	Medium	10" to 11" (25 to 28cm)
NC15433-3	Large	11" to 12" (28 to 30cm)
NC15433-4	X-Large	12" to 14" (30 to 35cm)



## NORCO™ FLOTATION ELBOW PAD

*Visco-elastic foam distributes pressure evenly.*

- Visco-elastic foam cushion with a knitted cotton, Orlon® and Lycra® sleeve.
- Use to protect the elbow from pressure sores or to lessen nerve hypersensitivity.
- Snug-fitting sleeve reduces migration of the pad on the arm.
- Pad measures 8" x 5½" x ½" (20 x 14 x 1.3cm) and sleeve measures 13" (33cm).
- Hand wash and air dry. Remove cushion prior to washing.
- Sold in pairs.
- Latex free.
- To size, measure circumference at elbow crease.

		Elbow Circumference
NC44889	Regular	9½" to 10½" (24 to 27cm)
NC44890	Large	10½" to 12" (27 to 30cm)



## SILIPOS® HEEL/ELBOW SLIPOVER

*Gel pad helps insulate against shock and vibration.*

- Tubular polyester sleeve with a polymer gel disc impregnated in the fabric.
- Use on the heel or elbow for comfortable pressure relief.
- Sold individually.
- Latex free.
- The small size measures 8" (20cm) long. The large size measures 10" (25cm) long.
- To size, measure circumference at elbow crease.

		Elbow Circumference
NC83022	Sm./Med.	11" to 13" (28 to 33cm)
NC83023	Lg./X-Lg.	13" to 15" (33 to 38cm)