



MEDI-RIP® TAPE

Self-adherent bandage provides comfortable breathable compression.

- Comfortable cotton tape offers excellent breathability and absorbency.
- Lightweight bandage adheres to itself without sticking to hair, skin or clothing.
- This non-constricting tape will stretch to almost twice its original length for excellent conformability on hard-to-bandage areas.
- Available in three widths. Rolls measure 5 yds. (4.6m) long.

	Width	
NC12601	1" (2.5cm)	24 Rolls
NC12602	2" (5.1cm)	12 Rolls
NC12603	3" (7.6cm)	12 Rolls



FABRIFOAM® NuStimWRAP™

Two-way stretch wrap offers non-slip support and compression.

- Patented fabric and foam composite material (not neoprene) that provides moderate elasticity.
- Use to wrap secure hot/cold packs and electrodes. Excellent for joint compression wraps and straps.
- Latex free.

Pkg. of 3

NC20092-1	2½" x 18" (6.4 x 46cm)
NC20092-2	2½" x 36" (6.4 x 91cm)
NC20092-3	2½" x 48" (6.4 x 122cm)
NC20092-4	4" x 18" (10 x 46cm)
NC20092-5	4" x 36" (10 x 91cm)
NC20092-6	4" x 48" (10 x 122cm)
NC20092-7	6" x 18" (15 x 46cm)
NC20092-8	6" x 36" (15 x 91cm)
NC20092-9	6" x 48" (15 x 122cm)



MOLINPIC TAPE

Co-adhesive tape provides excellent soft-tissue support.

- Self-adherent, elastic tape will not stick to skin, clothes or hair.
- Thicker, more dense than other cotton, self-adherent tapes.
- Ideal for use as flexible support when weaning off a rigid splint, recovering from sports injuries or when support for weak tendons is required.
- Available in red or blue. Each roll measures 2" x 4 yds. (5.1cm x 3.7m).

4 yd. (4.1m) Roll

NC12648-2	Red	2" (5.1cm)
NC12649-2	Blue	2" (5.1cm)



TUBIPAD® STOCKINETTE

Padding prevents skin breakdown.

- Elastic stockinette laminated to a foam lining.
- Helps control edema and prevent pressure sores.
- Provides maximum comfort while splinting clients who are extremely sensitive to heat.
- Sold in dispenser packs of 13 yd. (12m) rolls.
- Sizes measure flat width of stockinette.

	Flat Width	
NC15791	Small	3" (7.6cm)
NC15792	Medium	3½" (8.9cm)
NC15793	Large	4" (10cm)



ECONOMY COTTON ELASTIC BANDAGES

Economical, latex free bandages are a great value.

- These woven cotton bandages can be used for support and compression, or to secure dressings, ice, splints or casts.
- Wash with a mild detergent and warm water on a gentle cycle. Holds up with repeated washing.
- Each individually sealed bandage includes two elastic clips. Stretched length measures 4½ yds. (4.1m).
- Latex free.

Box of 10	Width
NC20083	3" (7.6cm)
NC20084	4" (10cm)
NC20086	6" (15cm)



PROFESSIONAL COTTON ELASTIC BANDAGES

Latex-free bandages provide controlled pressure.

- Similar to the Economy Cotton Elastic Bandages sold left, these Professional Bandages have about 30% less stretch.
- Woven cotton construction ensures improved comfort, reusability and breathability.
- Each individually sealed bandage includes two elastic clips.
- Stretched length measures 5 yds. (4.6m).
- Latex free.

Box of 10	Width
NC20072	2" (5.1cm)
NC20073	3" (7.6cm)
NC20074	4" (10cm)
NC20076	6" (15cm)

EDEMA CONTROL AND TAPING

DEMA® WRAP



Economical and lightweight, cohesive wrap provides excellent support and compression.

- Use to support sprains/strains and to apply compression for edema restriction.
- The selection of widths available ensure the ideal fit. Use the 1" (2.5cm) rolls for wrapping fingers, 2" and 3" (5.1 and 7.6cm) rolls best fit palms and wrists, while 4" (10cm) rolls are great for wrapping biceps and triceps.
- Use as a bandage and as a replacement for tape when securing dressings/catheters.
- Rolls measure 5 yds. (4.6m) long.

Standard Beige Width

NC68500	1" (2.5cm)	30 Rolls
NC68501	2" (5.1cm)	36 Rolls
NC68502	3" (7.6cm)	24 Rolls
NC68503	4" (10cm)	18 Rolls

Latex Free Beige Width

NC68510	1" (2.5cm)	30 Rolls
NC68511	2" (5.1cm)	36 Rolls
NC68514	3" (7.6cm)	24 rolls

Latex Free Color Pack Width

NC68512	1" (2.5cm)	30 Rolls
NC68513	2" (5.1cm)	36 Rolls

Latex Free Color Packs.



Make quick and easy light compression wraps.



No scissors required.



Provides compression and support as a soft restriction splint.



Self adheres and will not stick to skin, hair or clothing.

Secure splints while providing edema compression.

Adheres to itself, not to skin.



Use to provide circumferential support anywhere.



COBAN® SELF-ADHERENT WRAP

Wrap provides constant therapeutic pressure to help reduce edema.

- Coban® self-adherent elastic wrap is ideal for reducing edema and for holding splints, gel sheeting or elastomer molds in place.
- Reusable wrap provides light to firm compression.
- Wrap in a spiral fashion or figure 8 pattern to obtain desired compression where needed.
- Rolls are 5 yds. (4.6m) long.
- Latex-Free Coban® is available in beige rolls.

Beige	Width	
NC12651	1" (2.5cm)	30 Rolls
NC12652	2" (5.1cm)	36 Rolls
NC12653	3" (7.6cm)	24 Rolls
NC12654	4" (10cm)	18 Rolls
Blue	Width	
NC12646	1" (2.5cm)	30 Rolls
NC12647	2" (5.1cm)	36 Rolls

Latex-Free Beige	Coban® Width	
NC18215	1" (2.5cm)	30 Rolls
NC18216	2" (5.1cm)	36 Rolls



Coban® Self-Adherent Wrap can be used to stabilize elastomer splints.



Use Anchor™ Fix under Anchor™ Rigid Strapping Tape.

Anchor™ Rigid Taping Kit

ANCHOR™ RIGID STRAPPING TAPE AND ANCHOR™ FIX

Designed for corrective and preventative taping for shoulder, patellofemoral and lower extremity.

- **Anchor™ Rigid Strapping Tape** is a strong-adhering, rayon-backed, zinc oxide, brown fabric tape that can be used for correcting patella glide and other realignment techniques.
- **Anchor™ Fix** is a non-woven, hypoallergenic, white fabric underwrap that is used under Anchor™ Rigid Strapping Tape. Latex-free, acrylic adhesive helps to alleviate skin trauma when regularly taping with Anchor™ Rigid Strapping Tape.
- **Anchor™ Rigid Strapping Kit** includes one roll of 1½" (3.8cm) Anchor™ Rigid Strapping Tape, one roll of 2" (5.1cm) Anchor™ Fix, 10 Adhesive Remover Wipes (NC20418) and 10 Skin Prep/Barrier Wipes (NC20421).

Anchor™ Rigid Strapping Tape 15 yd. (14m) Rolls

NC22000-01	1½" (3.8cm)	1 Roll
NC22000-05	1½" (3.8cm)	5 Rolls
NC22000-30	1½" (3.8cm)	30 Rolls

Anchor™ Fix 10 yd. (9.1m) Rolls

NC22001-2	2" (5.1cm)	1 Roll
NC22001-2C	2" (5.1cm)	12 Rolls

Anchor™ Rigid Taping Kit

NC22003	(1)
---------	-----



SKIN PREP TAPE

Protects skin from tape irritation.

- Strong, waterproof and washable foam underwrap conforms to the body and accommodates tissue expansion.
- Each roll measures 2¾" x 30 yds. (7.0cm x 27m).
- Non-adhesive and latex free.

NC20700	Box of 12
NC20700-48	Box of 48



ATHLETIC TAPE

Offers great holding power for extremity taping.

- Made of strong, 100% bleached cotton fabric. Zinc oxide adhesive helps prevent irritation.
- Specially treated cloth backing helps prevent excessive perspiration entrapment.
- Non-elastic, rigid tape adheres strongly, but unwinds easily and consistently all the way to the core.
- Use for fingers, hands, wrists, elbows, knees and ankles.

15 yd. (14m) Rolls

NC35140-10	1" (2.5cm)	48 rolls
NC35140-15	1½" (3.8cm)	32 rolls
NC35140-20	2" (5.1cm)	24 rolls



ADHESIVE REMOVER WIPES

Clean surfaces for redressing.

- Gently removes adhesive residue.
- Contains skin-softening aloe.
- Nonsterile, 2-ply.

NC20418 (100)

SKIN PREP/BARRIER WIPES

Helps protect skin from irritants and contaminants.

- Forms a protective barrier between skin and adhesive products.
- Reduces the chance of excoriation.
- Nonsterile, 2-ply.

NC20421 (50)



LEUKOTAPE® P SPORTTAPE

Use as a positioning tape.

- High tensile-strength tape with rayon backing can be used to facilitate retraining of muscles.
- Ideal for patellofemoral and shoulder taping.
- Rubber-based adhesive holds tape securely in place so that the therapist can reposition tissue.
- Tape measures 1½" (3.8cm) wide.
- Sold in 15 yd. (14m) rolls.

NC20196-1	1 roll
NC20196-5	Box of 5 rolls
NC20196	Box of 30 rolls

SpiderTech™ Kinesiology Tape*

Engineered to mimic the thickness, weight and elasticity of human skin.

Scan to view
Elbow Spider video
and instructions



NEW SPIDERTECH™ TAPE

High-quality kinesiology tape provides superior adhesion with minimal irritation.

- Made of high-grade cotton material with 100% acrylic, hypoallergenic adhesive.
- “Gentle wave pattern” material provides superior adhesion with minimal irritation.
- Tape is water-resistant and breathable, and can be worn for up to five days.
- Latex free.

• **Rolls** measure 2" x 16.4' (5.0cm x 5m). Available in Blue, Black, Pink and Beige and **Gentle** (for sensitive skin).

• One **Bulk Roll** is equivalent to six standard rolls and is an economical solution for athletic or healthcare providers who use kinesiology tape consistently. Bulk Rolls makes taping convenient and help prevent the need to stop mid-application to retrieve another roll to finish taping. Bulk Rolls measure 2" x 103' (5.0cm x 31.5m).



Gentle tape is perfect for use with geriatric, pediatric and post-operative patients.



Rolls 2" x 16.4' (5.0cm x 5m)

NC22022-502	Beige	(6)
NC22022-503	Blue	(6)
NC22022-504	Red	(6)
NC22022-505	Black	(6)
NC22022-501	Gentle	(6)

Bulk Rolls 2" x 103' (5.0cm x 31.5m)

NC22021-502	Beige	(1)
NC22021-503	Blue	(1)
NC22021-504	Red	(1)
NC22021-505	Black	(1)
NC22021-501	Gentle	(1)

NEW SPIDERTECH™ X, Y, I POWERSTRIPS

Precut strips eliminate the need to cut the most commonly used kinesiology tape forms from rolls.

- X, Y or I Spiders can be applied independently or used in conjunction with SpiderTech™ pre-cut products to further enhance the therapeutic impact of treatments.
- X Y & I shapes replicate the most commonly used kinesiology tape therapy forms, providing the perfect “kit” to assemble complete kinesiology tape applications.
- “By the numbers” instructions on the backing simplify application and help to eliminate glue contamination.
- SpiderTech™ easy-tear perforated backing allows ease of use.
- Ideal as a precut solution for pediatric patients.
- Machine cut edges and radius curve ends enhance wearability and application endurance.
- Each box of PowerStrips is equal to two standard rolls of tape.

X and Y PowerStrips are sold in boxes of 40. **I PowerStrips** are sold in boxes of 100.

SpiderTech™ “I” PowerStrips

NC22047-501	Gentle	(100)
NC22047-503	Blue	(100)
NC22047-504	Red	(100)

SpiderTech™ “X” PowerStrips

NC22048-501	Gentle	(40)
NC22048-503	Blue	(40)
NC22048-504	Red	(40)

SpiderTech™ “Y” PowerStrips

NC22049-501	Gentle	(40)
NC22049-503	Blue	(40)
NC22049-504	Red	(40)



SpiderTech™
X, Y and I
PowerStrip
precut shapes.



SpiderTech™ precut kinesiology tape supports muscles for improved stability, security and power.

Use to improve athletic performance and for rehabilitation following surgery or injury. Helps restore normal muscle activation. Improves feedback and timing of muscles supporting the affected area. Relieves pressure at the points of pain while allowing the wearer to move freely. Provides protection during healing and aids in the prevention of injury. Each product has been


designed for use on pre-determined areas of the body, eliminating the need to design and cut the tape before applying. Can be applied in a multitude of ways with the 16 unique static designs that are currently available. Application efficiency saves valuable time for home or clinic use.

NEW SPIDERTECH™ ELBOW SPIDER

- Ideal for treating tennis elbow, golfer's elbow, and muscle strain.
- Sold in packs of 5.
- Latex free.

SpiderTech™ Elbow Spider

NC22027-502	Beige	(5)
NC22027-503	Blue	(5)
NC22027-504	Red	(5)
NC22027-505	Black	(5)
NC22027-501	Gentle	(5)



NEW SPIDERTECH™ WRIST SPIDER

- Ideal for treating sprained wrist, wrist instability, carpal tunnel syndrome, hypermobility and tendinopathies of the fingers.
- Sold in packs of 10.
- Latex free.

SpiderTech™ Wrist Spider

NC22028-502	Beige	(10)
NC22028-503	Blue	(10)
NC22028-504	Red	(10)
NC22028-505	Black	(10)
NC22028-501	Gentle	(10)

Wrist Spiders worn
two on each arm.



Visit www.ncmedical.com to see our complete selection of SpiderTech™ kinesiology taping products and instructional videos.

Two Wrist Spiders can be applied dorsally and volarly to treat carpal tunnel syndrome symptoms.

* All SpiderTech™ products are manufactured using the original Kinesiology Tape from Japan by Nitto Denko.

KINESIO® TEX GOLD™ TAPE

High-quality elastic tape designed for maximum range of motion.

- This durable, latex-free tape is water-resistant and can be worn up to four days without irritation, even during strenuous exercise or hydrotherapy. The heat-activated adhesive will not leave any residue when removed.
- Kinesio® Tape has 140% elasticity, the same basic elasticity found in skin and muscles. This allows the tape to work with the body, providing full support and stability without adversely affecting the natural range of motion.
- Ideal for patellofemoral treatment, shin splints, plantar fasciitis, rotator cuff injuries and carpal tunnel syndrome.

Kinesio® Tex Gold™ 16.4' (5m) rolls

NC89035-102	Beige	1" (2.5 cm)	(12)
NC89035-202	Beige	2" (5.1 cm)	(6)
NC89035-203	Blue	2" (5.1 cm)	(6)
NC89035-204	Red	2" (5.1 cm)	(6)
NC89035-205	Black	2" (5.1 cm)	(6)
NC89035-302	Beige	3" (7.6 cm)	(4)

Kinesio® Tex Gold™ 103.3' (31.5m) Bulk rolls

NC89035-332	Beige	2" (5.1 cm)	(1)
NC89035-333	Blue	2" (5.1 cm)	(1)
NC89035-334	Red	2" (5.1 cm)	(1)
NC89035-335	Black	2" (5.1 cm)	(1)



NEW KINESIO® FAN CUTS™

Minimize edema and inflammation with time-saving Kinesio® fan shaped precuts.

- **Fan Cuts™** are ideal for lymphatic draining and for treating edema and bruising. Useful for any condition that calls for increased circulation of blood and lymphatic fluid.
- Precut fan shape saves time. Rounded corners improve application wearability and longevity.
- 2" (5.1 cm) wide Fan Cuts™ feature five separate tape tails and are available in all four tape colors.
- 3" (7.6 cm) wide Fan Cuts™ feature eight separate tape tails and are only available in beige.
- Both sizes are cut to a 12" (30cm) length that may be cut down for specific patient taping.

Kinesio® Fan Cuts™, 2" (5.1 cm) width

NC89038-202	Beige	(50)
NC89038-203	Blue	(50)
NC89038-204	Red	(50)
NC89038-205	Black	(50)

Kinesio® Fan Cuts™, 3" (7.6 cm) width

NC89038-302	Beige	(50)
-------------	-------	------



KINESIO® PRO SCISSORS

Blades stay sharp for clean, easy cutting.

- High-carbon, stainless steel cutting surfaces remain sharp for years.
- Coated with black carbon and fluorine resin, for clean, quick and easy cuts without leaving residue.



NC89046

CLINICAL THERAPEUTIC APPLICATIONS OF THE KINESIO® TAPING METHOD

Corrective application techniques for advanced Kinesio® taping.

Describes six corrective application techniques in detail: mechanical, fascia, space, ligament/tendon, functional and lymphatic. Addresses specific treatment of clinical conditions. Soft cover, 198 pages with numerous photographs. By Kenzo Kase, DC, Jim Wallis, MS, ATC and Tsuyoshi Kase. Copyright 2003.

NC89047



CLINICAL KINESIO® TAPING DVD

Introduction, explanation and application of the Kinesio® Taping Method.

Certified Kinesio® Taping Instructor Jim Wallis, MS, ATC, demonstrates step-by-step application of Kinesio® Taping Techniques. Copyright 2000.

NC89048




KINESIO® TAPING FOR LYMPHOEDEMA AND CHRONIC SWELLING

Step-by-step instructions for over 108 applications.

By Kenzo Kase, DC, founder of the Kinesio® Taping Method, Kim Rock Stockheimer, MS, OTR/L, CKTI. Acknowledgement by Professor Neil Pillar. Copyright 2006.

NC60011



KINESIO® TAPING IN PEDIATRICS

Fundamentals and whole body taping for infants through adolescents.

By Kenzo Kase, DC, Founder of the Kinesio® Taping Method, Patricia Martin, PT, CKTI and Audrey Yasukawa, MOT, OTR/L, CKTI. Copyright 2006.

NC60012

