



NEW MEPILEX® FOAM DRESSING

Minimize both trauma to the wound and pain during dressing changes.

- Soft, conformable silicone foam dressing that absorbs exudate and maintains a moist wound environment.
- Safetac technology prevents dressing from sticking to the wound bed, yet allows gentle adherence to the surrounding dry skin.
- May be lifted and adjusted without losing its adherent properties.
- Minimizes the risk of maceration.
- Can be used on moderately exuding wounds, traumatic wounds, neuropathic ulcers, graft and donor sites.
- Sterile.

NC65448 4" x 4" (10 x 10cm) (5)



KERLIX® ROLLS

Unique soft gauze provides cushioning and conformability.

- Ideal for bandaging heads, limbs and difficult-to-dress wounds such as burns, plastic and orthopedic wounds.
- Made of pre-washed, fluff-dried, 100% woven gauze with unique crinkle-weave pattern provides loft and bulk to cushion and protect wound areas.
- Six-ply construction with finished edges reduces loose ends and lint.
- Sterile rolls are individually wrapped.
- Latex free.

2¼" x 3 yds. (5.7cm x 2.7m)

NC20131 (12)

NC20131C (96)

4½" x 4 yds. (11cm x 3.7m)

NC20132 (10)

NC20132C (100)



MEPISEL®

Helps make bandage changes less painful.

- Non-adherent silicone dressing in a fine mesh netting, coated with medical-grade silicone gel.
- Will not stick to moist tissue such as a wound bed, but adheres firmly to intact skin.
- Ideal for use with skin grafts, second degree burns and granulating wounds.
- Can be left over the wound bed for several days, while the outer dressing is changed as needed.
- Latex-free.

NC65503 2" x 3" (5 x 7.5cm) Box of 10

NC65504 3" x 4" (7.5 x 10cm) Box of 10

NC65508 4" x 7" (10 x 18cm) Box of 10



ELASTO-GEL™ OCCLUSIVE DRESSING

Soothes wounds and helps keep them clean.

- Sterile, opaque, absorbent hydro-gel pads.
- Helps reduce the pain and trauma of burns and wounds.
- Assists in keeping wounds free of contamination.
- Gel pad retains wound drainage preventing transfer to outer dressing layers. Latex-free.

NC20650 4" x 4" (10 x 10cm) Box of 25

NC20651 6" x 8" (15 x 20cm) Box of 5



N-TERFACE

Ideal for post-operative wound care.

- Lightweight, porous, second-skin material prevents dressing from adhering to wounds. Light tack keeps the material in position but will not pull on skin or wounds when removed.
- Can be repositioned as needed.
- Sterile, individually wrapped dressings Latex-free.

NC20009 4" x 4" (10 x 10cm) (25)

NC20010 4" x 12" (10 x 30cm) (12)



SPONGES

Versatile sponges provide maximum absorption.

- Cotton-filled sponges can be used for a wide range of dressing, cleaning and prepping procedures.
- Cut edges are folded—no rough edges are exposed.
- Non-sterile Sponges are available in two styles: All-gauze and Non-woven.
- Sterile sponges individually wrapped.
- Latex-free.

All-gauze, non-sterile sponges: Box of 200

NC20270-208 2" x 2" (5.1 x 5.1cm) 8-ply

NC20270-308 3" x 3" (7.6 x 7.6cm) 8-ply

NC20270-408 4" x 4" (10 x 10cm) 8-ply

NC20270-412 4" x 4" (10 x 10cm) 12-ply

Non-woven, non-sterile sponges: Box of 200

NC20273-204 2" x 2" (5.1 x 5.1cm) 4-ply

NC20273-304 3" x 3" (7.6 x 7.6cm) 4-ply

NC20273-404 4" x 4" (10 x 10cm) 4-ply

All-gauze, sterile sponges: Box of 25

NC20166 3" x 3" (7.6 x 7.6cm) 12-ply

NC20167 4" x 4" (10 x 10cm) 12-ply



GAUZE SPONGES

Open weave maximizes conformity.

- Sterile, open-weave, 12-ply dressings are made from high-grade cotton and measure 4" x 4" (10 x 10cm). Exceeds USP Type VII gauze requirements.
- Use for wound dressing, general cleaning, minor prepping, wound packing and debriding.
- Latex free.

Sterile (Individually wrapped, 12-ply sponges)

NC20163K 4" x 4" (10 x 10cm) (50)

WOUND CARE/DRESSINGS AND BANDAGES



XEROFORM® PETROLATUM DRESSING

Clings and conforms to body contours.

- Non-adherent, 3% Bismuth Tribromophenate in a non-medicinal, petrolatum blend on fine mesh gauze. Bacteriostatic action.
- Prevents gauze from adhering to wounds.
- Individually wrapped in sterile packaging.
- Latex free.

Box of 50

NC20141 1" x 8" (2.5 x 20cm)

NC20142 5" x 9" (13 x 23cm)



CURASORB™ CALCIUM ALGINATE®

Absorbs approximately 20 times the weight of the dressing.

- Absorbs exudate from the wound and creates a protective cushioning gel and optimal moist healing environment. Protective gel binds exudate and helps prevent skin maceration.
- Easy, painless dressing removal minimizes tissue trauma.
- Individually wrapped in sterile packaging.
- Latex free.

Package of 5

NC95206-1 12" (30cm) rope

Package of 10

NC95206-2 2" x 2" (5.1 x 5.1cm)

NC95206-3 4" x 4" (10 x 10cm)



MEFIX®

Use under sports tape or to secure dressings.

- Non-woven, 100% polyester fabric with a non-irritating adhesive layer. Air permeable.
- Porous and non-absorbing. Minimizes fluid strike-through. Use to affix gauze, tubes, cannulas, topical dressings and drainage devices.
- Stays in place, yet can be removed gently.
- Latex free.

11 yd. (10m) Rolls

NC65210 1" (2.5cm)

NC65220 2" (5.1cm)

NC65240 4" (10cm)

NC65260 6" (15cm)



ULTEC PRO®

ALGINATE HYDROCOLLOID DRESSINGS

Flexible dressing for dry to nearly dry wounds.

- Highly conformable dressing designed for use on lightly exuding to dry wounds.
- Creates and maintains an optimal moist healing environment.
- Unique formulation absorbs exudate and helps prevent the exudate from causing maceration of the peri-wound skin.
- Individually wrapped in sterile packaging.
- Latex free.

Regular (Package of 5)

NC95204-1 4" x 4" (10 x 10cm)

Sacral (Package of 5)

NC95204-4 4" x 5" (10 x 13cm)



TENDERWRAP® UNNA BOOT BANDAGE

Non-hardening, zinc oxide formulation.

- Soft, non-raveling, 100% cotton gauze is impregnated with a non-hardening zinc oxide paste.
- Use for venous stasis conditions.
- Bandage remains flexible, encouraging patient ambulation.
- Conforms well to body contours and provides even compression.
- Gelatin-free, preservative-free.
- Individually wrapped in non-sterile packaging.
- Latex free.

3" x 10 yds. (7.6cm x 9.1m)

NC95209 (1)

NC95209-3 (12)

4" x 10 yds. (10cm x 9.1m)

NC95214 (1)

NC95209-4 (12)

TELFA®

ISLAND DRESSINGS

Convenient, all-in-one dressing.

- Ideal dressing for use on lightly draining wounds.
- Soft, nonwoven backing conforms around wound and seals on all four sides.
- Provides protection from external environment. Non-adherent TELFA® will not disrupt healing tissue.
- A secure dressing for central and peripheral IV sites.
- Individually wrapped dressings in a sterile, peel-back package.
- Latex-free.

Package of 25

NC95202-2 4" x 5" (10 x 13cm)

NC95202-3 4" x 8" (10 x 20cm)



CURITY®

NON-ADHERING DRESSINGS

Allows drainage of exudate away from wound.

- Oil-emulsion blend impregnated into an open-mesh, knitted fabric dressing.
- Ideal for donor sites, delicate suture lines, graft sites and draining wounds.
- Individually wrapped in sterile packaging.
- Latex-free.

Package of 50

NC95213-2 3" x 3" (7.6cm x 7.6cm)

Package of 24

NC95213-3 3" x 8" (7.6cm x 20cm)





ANCHOR™ FIX

Use this underwrap adhesive tape to secure bandages on sensitive skin.

- Anchor™ Fix is a non-woven, hypoallergenic, white fabric tape. Conforms well to bony areas.
- Acrylic adhesive is better tolerated by skin and helps to alleviate skin trauma.
- Latex free.



10 yd. (9.1m) Rolls

NC22001-2	2" (5.1cm)	1 Roll
NC22001-2C	2" (5.1cm)	Box of 12



CO-WRAP™ SELF-ADHERENT TAPE

A variable alternative to bandages or tape.

- Made of 90% cotton and 10% polymer.
- Slight tackiness helps keep it adhered to itself.

5 yd. (4.6m) Rolls

NC12611	1" (2.5cm)	24 rolls
NC12612	2" (5.1cm)	12 rolls
NC12613	3" (7.6cm)	12 rolls



TRANSPORE™ ADHESIVE TAPE

Adheres securely, even when wet.

- Transparent, plastic-like, surgical tape.
- Strong adhesive effectively holds electrodes on the skin.
- Waterproof.
- Latex free.

Box of 12

NC20032	1" x 10 yds. (2.5cm x 9.1m)
---------	-----------------------------



COVER-ROLL® STRETCH

Covers an entire wound with one application.

- Non-woven, air permeable, self-adhesive and hypoallergenic. Latex free.
- Stretches and conforms to bony areas.

2" x 2 yds. (5.1cm x 1.8m)

NC20192	(1)
NC20192-12	(12)

2" x 10 yds. (5.1cm x 9.1m)

NC20194	(1)
NC20194-12	(12)

4" x 2 yds. (10cm x 1.8m)

NC20190	(1)
NC20190-12	(12)

4" x 10 yds. (10cm x 9.1m)

NC20191	(1)
NC20191-12	(12)



CURI-STRIP® ADHESIVE WOUND CLOSURES

Hypoallergenic adhesive helps keep wounds closed.

- Conformable, spunbonded nylon backing free of reinforcing filaments for increased freedom of movement.
- Helps minimize premature pop-off and formation of stress blisters.
- Highly porous backing for increased oxygen exposure. Superior moisture vapor transmission for minimal skin maceration.
- Latex free.

5 per envelope, 50 envelopes per box

NC95200-1	1/8" x 3" (3.2mm x 7.6cm)	(250)
-----------	---------------------------	-------

6 per envelope, 50 envelopes per box

NC95200-2	1/4" x 1 1/2" (6.4mm x 3.8cm)	(300)
NC95200-5	1/2" x 4" (13mm x 10cm)	(300)



MICROFOAM™ TAPE

Thick, padded tape with light adhesive is great for wound care dressing.

- Closed-cell padding in a pressure-sensitive tape.
- Light adhesive allows easy removal without residue.
- Use as a liner for small hand or finger splints.
- May be used as a compressive wrap to treat edematous fingers.
- Cushions edges of splints or equipment.
- Water resistant.
- Latex free.

Package of 12 Rolls

NC12640	1" x 5 yds. (2.5cm x 4.6m)
---------	----------------------------

MICROPORE® TAPE

Hypoallergenic tape is great for sensitive skin.

- Thin, lightweight cotton-paper surgical tape is easy to cut or tear by hand. Ideal for edema control with hypersensitive skin.
- Tiny holes throughout increase breathability and aid in the healing process.
- Low-tack adhesive will not irritate skin during removal. A great choice for frequent dressing changes. Latex free.



Box of 12 Rolls

NC20007-1	1" x 10 yds (2.5cm x 9.1m)
-----------	----------------------------

Box of 6 Rolls

NC20007-20	2" x 10 yds (5.1cm x 9.1m)
------------	----------------------------

WOUND CARE/BANDAGES



WEBRIL® UNDERCAST PADDING

Ideal for holding dressings in place or for use under a cast or splint.

- Made of non-sterile cotton material.
- Latex free.

4 yd. (3.7m) Rolls

NC20027	2" (5.1cm)	24 rolls
NC20028	3" (7.6cm)	12 rolls
NC20029	4" (10cm)	12 rolls



COVERLET®

Sterile, non-constricting dressing with adhesive edges conforms to fingertips and knuckles.

- Allows free movement at the DIP joint.
- Latex free.

NC20021	1" x 3" (2.5 x 7.6cm)	(100)
NC20022	Knuckle	(100)
NC20023	Small Digit	(100)
NC20024	Large Digit	(50)



SURGITUBE®

Keeps dressings securely on fingers.

- This tubular gauze is made of 100% cotton.
- Provides maximum aeration. Secures dressings.
- An applicator is included in the 5 yd. (4.6m) box only.
- Latex free.

NC15797	3/8" x 5 yds. (16mm x 4.6m)
NC15798	7/8" x 5 yds. (22mm x 4.6m)
NC15797-50	3/8" x 50 yds. (16mm x 46m)
NC15798-50	7/8" x 50 yds. (22mm x 46m)



CONFORM™ STRETCH BANDAGES

Provides softness, conformability, low lint and high absorbency.

- One-ply, cotton/polyester blend crocheted bandage. Holds securely to any body contour; allows for movement and some soft tissue swelling.
- Sterile or non-sterile.
- Ideal for securing dressings, IVs and splints or for providing mild compression or support. Stays in place with minimal taping.
- Rolls measure 4.1 yds. (3.7m) stretched.
- Latex free.

Sterile

NC95208-1	1" (2.5cm)	12 Rolls
NC95208-2	2" (5.1cm)	12 Rolls
NC95208-3	3" (7.6cm)	12 Rolls

Non-sterile

NC95208-6	1" (2.5cm)	24 Rolls
NC95208-7	2" (5.1cm)	12 Rolls
NC95208-8	3" (7.6cm)	12 Rolls

DISPOSABLE SCALPELS

Sterile scalpels have a stainless steel blade and a plastic handle.

- Individually wrapped.
- Designed for one-time use.

Package of 10
NC20160-10



NEW PEHA-HAFT® LF

Cohesive conforming gauze bandage.

- Ideal for use as a first aid dressing, and in applications when light and consistent compression is required. Can be used as a secondary dressing or retention bandage.
- Offers a combination of light absorption, gentle support and conformability.
- Reduces risk of sensitivity, especially on delicate or traumatized skin.
- Adheres to itself, never to skin or hair, without pins, clips or tape.
- Latex free.

1" (2.5cm) Width

NC20225-01	4.5 yds. (4.1m)	(8)
------------	-----------------	-----



TUBITON® APPLICATOR

For easy application of tubular bandages.

- Can be used with Tubegaуз® Finger Bandage and BandNet™ (shown above right).
- Stretches and holds the bandage open for comfortable donning. Minimizes skin contact by preventing the bandage from rubbing against the finger when applied.

NC15755



BANDNET™

Tubular netting secures bandages without tape.

- Made of expandable nylon webbing.
- Each box contains 25 yds. (23m).
- Use the Tubiton® Applicator® (shown below left) for easy application. Sold separately.

Size Fits

NC20013	1	Fingers, Toes
NC20014	2	Lg. Fingers, Wrists
NC20015	3	Knees, Feet, Elbows, Hands
NC20016	4	Knees, Feet, Elbows, Hands
NC20017	5	Knees, Feet, Elbows, Hands



TUBEGAUZ® FINGER BANDAGE

Ideal for holding dressings on fingers.

- Made of 100% white seamless gauze cotton.
- Use the Tubiton® Applicator (shown left) for easy application. Sold separately.

NC20019	3/8" x 50 yds. (16mm x 46m)
---------	-----------------------------



SUTURE REMOVAL KIT

Disposable kits ensure sterility.

- Sterile kit contains metal-insert forceps, Littauer Scissors and a cotton-filled gauze sponge.
- Ideal for wound care following injury or surgery.
- Designed for one-time use.

NC20149 Single kit
NC20150 Box of 50 kits



COTTON-TIPPED APPLICATORS

Soft and absorbent tips.

- Single-tipped wood stick applicators.
- Use for wound debridement or application of topical agents.
- Non-sterile applicators measure 6" (15cm) long and are bulk-packaged.
- Sterile applicators measure 6" (15cm) long and are individually wrapped.

NC95215-6 Non-sterile (100)
NC20050 Sterile (100)

HIBICLENS®

Reduce the risk of cross infection.

- Hibiclens® can be used as a skin cleanser in wound care programs.
- For external use only. Do not use near eyes, or mouth.

4 fl. oz. (118ml)

NC12690

32 fl. oz. (946ml)

NC12690-1



EXAM GLOVES

Protect hands during treatment and therapy.

- Sterile gloves are individually wrapped and sold in pairs.

Nitrile Non-Sterile, Latex Free and Powder Free: Box of 100

NC20510-1 Small
NC20510-2 Medium
NC20510-3 Large
NC20510-4 X-Large

Vinyl Non-Sterile, Latex Free: Box of 100

NC20517-1 Small
NC20517-2 Medium
NC20517-3 Large
NC20517-4 X-Large

Vinyl Sterile, Latex Free: Box of 50

NC20572-1 Small
NC20572-2 Medium
NC20572-3 Large

Vinyl Non-Sterile, Powder Free: Box of 100

NC20520-1 Small
NC20520-2 Medium
NC20520-3 Large
NC20520-4 X-Large

Latex Sterile: Box of 50

NC20511-0 Small
NC20511-1 Medium
NC20511-2 Large

Latex Non-Sterile: Box of 100

NC20518-1 Small
NC20518-2 Medium
NC20518-3 Large
NC20518-4 X-Large

Latex Non-Sterile and Powder Free: Box of 100

NC20519-1 Small
NC20519-2 Medium
NC20519-3 Large
NC20519-4 X-Large



BANDAGE SCISSORS

Perfect for trimming fracture braces.

- Stainless steel blades measure 2½" (6.4cm).
- Scissors feature blunt safety tips.
- Autoclavable.
- For right-handed use only.

NC12715



LISTER BANDAGE SCISSORS

Two sizes to suit different needs.

- Blunt-sided utility scissors are made of stain-resistant steel to resist corrosion and rusting.
- Ideal for removing gauze or bandages
- Blunt end decreases the chance of accidental skin injuries.
- Autoclavable.
- For right-handed use only.

	Total Length
NC12742	5½" (14cm)
NC12742-1	7¼" (18cm)



SKIN PREP/BARRIER WIPES AND ADHESIVE REMOVER WIPES

Skin Prep/Barrier Wipes protect skin; Adhesive Remover Wipes prepare skin for redressing.

- **Skin Prep/Barrier Wipes** form a protective barrier between skin and adhesive products. Helps reduce the chance of excoriation. Protects skin from feces, urine, chemical irritants and gastrointestinal fluids. Large, two-ply, non-sterile wipes.
- **Adhesive Remover Wipes** provide a clean surface for redressings. Contains aloe to leave the skin soft after usage.

NC20421 Skin Prep/Barrier Wipes (50)
NC20418 Adhesive Remover Wipes (100)