



— Air glide roller ball glides smoothly over bare skin or clothing, eliminating the use of messy gels and oils.

— Ergonomically designed handle reduces fatigue.

Model 610

Model 300



— Acupressure endpoint delivers precise compression to sore trigger points.

GLIDEPOINT MASSAGE TOOLS

Ergonomically designed tools can be applied to soft tissue for longer periods of time without creating itching or numbness.

The air glide roller ball on the **Glidepoint 300** moves smoothly over bare skin or clothing without requiring massage gels and oils. This model is ideal for foot self-massage. When the Glidepoint 300 is placed on the floor, the user's body weight applies pressure with a minimum of effort.

Glidepoint Wand Massager 610 has an extended, textured handle for easy grip and maneuverability. Ergonomically designed handle reduces hand fatigue and facilitates greater massage pressure. An acupressure endpoint on the handle is ideal for precise compression of sore trigger points.

NC92208-1 Model 300

NC92208 Model 610



BACKNOBBER® II

Simple, deep-muscle self-massage tool.

This deep-muscle massage tool helps people manage chronic muscular pain on their own or in conjunction with professional medical treatment and therapy. The Backnobber® II is designed to be hooked over the shoulder or under the arm so that one of the two ends can be pressed into the muscles of the upper or lower back or neck. The Backnobber® II's natural shape and dimension allow the user to apply as much pain-releasing pressure for as long as required. Made of a highly durable, polymer composite molded in the shape of an "S". A locking mechanism allows it to be taken apart for convenient storage or travel. Fits into an average size briefcase or overnight bag.

NC92202



Knobs serve as hand holds for self-applied massage.


THERA CANE®

Ideal tool for self-management of soft tissue syndromes.

Use Thera Cane® to treat trigger points associated with myofascial pain syndrome. The curved, fiberglass cane has several knobs for easy use. Little force is required to apply pressure to trigger points. An instructional video is available separately.

NC92220 Thera Cane®

NC92221 Video



JACKNOBBER® II

Applies deep pressure.

Help relieve muscle stiffness, fatigue and painful spasms with the Jacknobber II. Two small knobs and two large knobs offer precise control. Place on the floor or against a wall and lean on the upright knob for different massaging exercises. Instructions included. Color may vary.

NC92201

NEW! OAKWORKS® PORTABLE MASSAGE TABLES

Available in One™ (economy) and Aurora™ (standard) models to suit your needs.

These durable, lightweight massage tables feature:

- Patented cable lock system constructed of quality hardwood and aircraft tested cables.
- Wooden Easy-Adjustment Knobs with oil- and dirt-resistant grips.
- High-density, multi-layer semi-firm padding for longevity and comfort.
- Super soft, earth-friendly TerraTouch™ upholstery in Sage or Ocean Blue color as indicated.
- Durable, heavy-duty leg extensions.
- Call for freight quote.

Oakworks® One™ Table

Available in Ocean Blue (shown) and Sage.



Aurora™ table

Available in Sage.

Classic Face Rest
Available in Sage.



NEW! OAKWORKS® ONE™

Economy model is one incredible value.

Portable table features face rest outlets at both ends. Includes arm hammock, adjustable face rest and pad, access end panel, shiatsu cables and carrying case with shoulder strap and large zippered pocket.

- Length x Width: 73" x 30" (185 x 76 cm).
- Height: 24" to 34" (61 to 86 cm).
- Weight: Approximately 34 lbs. (15 kg).
- CA 117 compliant.
- Weight capacity: 500 lbs. (227 kg).

NC82091 Sage

NC82091-1 Ocean Blue

NEW! OAKWORKS® AURORA™

Durable, lightweight table transports easily.

- Length x Width: 73" x 30" (185 x 76 cm).
- Height: 24" to 34" (61 to 86 cm).
- Weight: 25 to 29 lbs. (11 to 13 kg).
- Weight capacity: 550 lbs. (249 kg).

NC82090 Sage

NEW! CLASSIC FACE REST

Adds up to 12" (30 cm) to portable tables.

Provides excellent access to head, neck and shoulder while prone or supine. Strong, lightweight aluminum. Durable, powder-coated scratch-resistant finish. Includes Memory Plus Crescent face rest.

NC82093 Sage



BIOTONE® DUAL-PURPOSE MASSAGE CREAM

Combines the benefits of both oil and lotion.

Dual-purpose, silky massage cream has a rich, smooth consistency for even application. Remains effective without frequent re-application. Softens and moisturizes skin without a greasy feel.

Tube

NC70234 7 fl. oz. (207 ml)

Jars

NC70231 14 fl. oz. (414 ml)

NC70232 36 fl. oz. (1 liter)

NC70235 1 gallon (3.8 liter)



FREE-UP® SOFT TISSUE MASSAGE CREAM


Superb medium is formulated exclusively for healthcare professionals.

Excellent soft-tissue medium for surface or deep-tissue massage. Provides great glide, promoting exceptional tissue sense to fingers and hands. Safe, quality formula is bacteriostatic. Does not contain beeswax or fragrance. Not oily or slick. Does not inhibit deep tissue perception or discrimination.

NC70227 Unscented 16 fl. oz. (473 ml)

NC70228 Scented 8 fl. oz. (236 ml)

NC70229 Scented 16 fl. oz. (473 ml)



PRON PILLŌ™

Support the head, neck and shoulders during massage.

Increase comfort during massage by elevating and supporting the head, neck and shoulders in a natural position. Durable pillow conforms to body contours yet returns to its original shape. Measures 19" long x 20" wide x 4" deep (48 x 51 x 10 cm). Cutout measures 3 1/2" wide x 8" long (8.9 x 20 cm) and the lip measures 3" (7.6 cm) wide. May be cleaned with disinfectant.

NC92343



PRELIM™ PROFESSIONAL GRADE BALM, LOTION AND CREAMS

Selection of formulas are ideal for band therapy and scar tissue massage.

The entire Prelim™ family of hypoallergenic massage products is free of fragrance and beeswax, the perfect choice for people with sensitive or hypersensitive skin. Each of these products has its own unique characteristics and properties. Choose the formula that best suits specific needs. Prelim™ is available in **Balm**, **Lotion**, **Soft Tissue Cream** and a **Dual-Action Cream** that is specially designed for both soft-tissue massage and deep-tissue mobilization.

PRELIM™ BALM

Creates surface friction for massaging large areas.

Contains carnauba, coconut oil and lanolin for moisturizing. Special formula creates surface friction so that large areas or muscle groups can be massaged without frequent re-application. Ideal for use with scar tissue massage tools such as probes, roller ball sticks or vibrators.



Balm

NC31201 15 oz. (425g)

PRELIM™ LOTION

Softens, moisturizes and replenishes skin.

This penetrating, non-greasy, professional-grade therapeutic lotion is formulated with cocoa butter, mineral oil, glycerin and vitamin E. Recommended for scar tissue and retrograde massage, and for use on burn, trauma and surgical scars. Squeeze tube allows for easy dispensing.



Lotion

NC31202 2 oz. (57g)

PRELIM™ SOFT TISSUE CREAM

Moisturize and massage scar tissue.

Prelim™ Soft Tissue Cream is ideal for surface or deep-tissue massage. Moisturizes the skin and provides glide and tissue perception over scar tissue and adhesions without friction, oiliness or slickness.

Soft Tissue Cream

NC31203 16 oz. (454g)



NEW! PRELIM™ DUAL-ACTION CREAM

Easy-glide formula massages and moisturizes.

Extra smooth formula applies like a rich cream and works like an oil, achieving the performance benefits of both textures. Provides easy glide and a consistency that offers excellent workability and absorbency. Softens and moisturizes the skin with no greasy feel.



Dual-Action Cream

NC31198 15 oz. (425g) jar

Heat attachment



WAHL® COMFORT VIBRATORS

Heat adds therapeutic comfort to massage.

Includes seven attachments to fit a variety of body contours. The Wahl® 8 Way Comfort Vibrator also includes a heating head attachment (replacement heat attachment sold separately). Two speeds: 60 cycles per second at low speed and 120 cycles per second at high speed. Both of these lightweight units weigh 1 lb. (.45kg). and measure 6½" (17cm) in length. Durable nylon casing houses a quiet, 120 volt, 16 watt motor. Fits in the hand comfortably.

NC70203 Wahl® 8 Way Vibrator w/ Heat

NC70203-1 Replacement Heat Attachment

NC70204 Wahl® 7 Way Vibrator



NEW! PROSSAGE™ HEAT


Provides gentle heat and area specific, controlled glide.

Prossage™ Heat is an all-natural, warming medium that consists of safflower seed oil, menthol, lanolin and lavender oil. Specially formulated for deep-tissue work, myofascial release and trigger point therapy. The combination of heat and friction increases blood flow and oxygen, reducing spasms. Prossage™ Heat makes it easier to hook the deep fascia, allowing the professional to work faster and more efficiently. Soft-tissue release is deeper and faster, resulting in muscle relaxation at the joints and tendons, reducing or eliminating trigger points, myofascial pain, tendinitis and cartilage degeneration due to progressive joint arthritis.

NC70067 3 oz. (85g)

NC70067-1 8 oz. (227g)

NC70067-2 32 oz. (907g)




BEARD'S MASSAGE

Guide to soft-tissue massage.

Topics include stroking, pressure and percussion manipulations, draping and positioning, equipment and more. Contains 320 descriptive illustrations with directional arrows. Soft cover, 185 pages. By Giovanni De Domenico and Elizabeth C. Wood. Fourth edition, copyright 1997.

NC35120



URIAS® AIR SPLINTS



Transparent air splints facilitate correct positioning and support compression.

Made of durable PVC sheeting with an inner layer that provides superior molding around the limb. Latex-free splints are open at both ends. Plastic zipper closure for easy donning and doffing. Use to immobilize the distal joints and to restore motion and control at the proximal joints. Apply to limb and inflate using a filter and a disposable mouthpiece (sold separately). Filters help prevent moisture from accumulating in splints. Filters change color to indicate the need for replacement. Clean with warm, soapy water. To size, measure the length of the extremity. Measurement equals the length of the splint before inflation.

Adult		Splint Length	Chambers
NC13004	Hand	8" (20 cm)	Double
NC13009	Hand/Wrist	12" (30 cm)	Double
NC13002	Half Arm	21" (53 cm)	Single
NC13001	Short Arm	27" (69 cm)	Single
NC13000	Long Arm	31" (79 cm)	Single
NC13003	Elbow	14" (36 cm)	Single
NC13013-60	Small Leg (no foot)	22" (56 cm)	Double
NC13013-70	Large Leg (no foot)	26" (66 cm)	Double
NC13013-F	Foot	12½" (32 cm)	Single
Child		Splint Length	Chambers
NC13011-20	Infant Arm	8" (20 cm)	Single
NC13011-30	Infant Arm	12" (30 cm)	Single
NC13011-40	Infant Arm	15¾" (40 cm)	Single
NC13012-20	Child Arm	8" (20 cm)	Single
NC13012-30	Child Arm	12" (30 cm)	Single
NC13012-40	Child Arm	15¾" (40 cm)	Single
NC12999	Child Arm	19¾" (50 cm)	Single
NC13012	Child Hand	8" (20 cm)	Double
NC13014-30	Child Leg (no foot)	12" (30 cm)	Double
NC13014-40	Child Leg (no foot)	15¾" (40 cm)	Double
NC13005	Mouthpieces (100)		
NC13006	Filters (5)		




FOAM ARM SUPPORT

Assist edema control and protect flaccid extremities.

This contoured, solid foam support is ideal for positioning the arm during treatment and exercise. Positions and elevates to 5" (13 cm) on a lapboard, table or bed without slipping. Support measures 20" (51 cm) long and 5" (13 cm) high at the tall end.

NC94133




VINYL COATED FOAM ARM SUPPORT

Support can be sterilized for clinic use.

Ideal for positioning arms during treatment and exercise. Vinyl coating resists stains and can be washed with anti-bacterial soap. Support measures 20" (51 cm) long and 5" (13 cm) high at the tall end.

NC94136





ANCHOR™ RIGID STRAPPING TAPE

Corrective and preventative taping for shoulder, patellofemoral and lower extremity.

Anchor™ Rigid Strapping Tape is a strong-adhering, rayon-backed, zinc oxide, brown fabric tape that can be used for realignment techniques like correcting patella glide. **Anchor™ Fix** is a non-woven, hypoallergenic, white fabric underwrap that is used under Anchor™ Rigid Strapping Tape when the tape must be secure under extreme stress and/or sweat, as with knee or shoulder taping. Latex-free, acrylic adhesive is better tolerated by skin and helps to alleviate skin trauma when regularly taping with Anchor™ Rigid Strapping Tape. **Anchor™ Rigid Strapping Kit** includes one roll of 1½" (3.8cm) Anchor™ Rigid Strapping Tape, one roll of 2" (5.1cm) Anchor™ Fix, 10 Adhesive Remover Wipes (NC20418) and 10 Skin Prep/Barrier Wipes (NC20421).

Anchor™ Rigid Strapping Tape



Anchor™ Rigid Strapping Tape 15 yd. (14m) Rolls

NC22000-01	1½" (3.8cm)	1 Roll
NC22000-05	1½" (3.8cm)	5 Rolls
NC22000-30	1½" (3.8cm)	30 Rolls

Anchor™ Fix 10 yd. (9.1m) Rolls

NC22001-2	2" (5.1cm)	1 Roll
NC22001-2C	2" (5.1cm)	12 Rolls
NC22001-4	4" (10cm)	1 Roll
NC22001-4C	4" (10cm)	12 Rolls

Anchor™ Rigid Taping Kit

NC22003	(1)
NC22003-12	(12)



TAPING TECHNIQUES: PRINCIPLES AND PRACTICE

Practical handbook of functional taping techniques.

Second edition covers taping techniques for treatment of muscle imbalances, unstable joints and neural control. New techniques include Mulligan and proprioceptive taping and other recent innovations. Also covers stabilization techniques and the theory and principles of unloading painful structures. Softcover, 240 pages. Copyright 2004. By Rose MacDonald.

NC20705



COVER-ROLL® STRETCH

Use as an underwrap tape or use to cover wounds.

Bandage is non-woven, air permeable, self-adhesive and hypoallergenic. Stretches and conforms to bony areas. Available in 2" or 4" (5.1 or 10cm) widths. Sold in 2 or 10 yd. (1.8 or 9.1m) rolls. Latex free.

2 yd. (1.8m) Rolls

NC20192	2" (5.1cm)	(1)
NC20192-12	2" (5.1cm)	(12)
NC20190	4" (10cm)	(1)
NC20190-12	4" (10cm)	(12)

New! 10 yd. (9.1m) Rolls

NC20194	2" (5.1cm)	(1)
NC20194-12	2" (5.1cm)	(12)
NC20191	4" (10cm)	(1)
NC20191-12	4" (10cm)	(12)



LEUKOTAPE® P SPORTSTAPE

Use as a positioning tape.

High tensile-strength tape with rayon backing can be used to facilitate retraining of muscles. Ideal for patellofemoral and shoulder taping. Rubber-based adhesive holds tape securely in place so that the therapist can reposition tissue. Tape measures 1½" (3.8cm) wide. Sold in 15 yd. (14m) rolls.

NC20196-1	1 roll
NC20196-5	Box of 5 rolls
NC20196	Box of 30 rolls



ATHLETIC TAPE

Great holding power for extremity taping.

Made of strong, 100% bleached cotton fabric. Zinc oxide adhesive helps prevent irritation. Latex free. Specially treated cloth backing helps prevent excessive perspiration entrapment. Non-elastic, rigid tape adheres strongly, but unwinds easily and consistently. Rolls measure 15 yds. (14m).

15 yd. (14m) Rolls

NC35140-10	1" (2.5cm)	48 rolls
NC35140-15	1½" (3.8cm)	32 rolls
NC35140-20	2" (5.1cm)	24 rolls

Skin Prep Wipes



Adhesive Remover Wipes



SKIN PREP/BARRIER WIPES

Helps protect skin from irritants and contaminants.

Forms a protective barrier between the skin and adhesive products. Reduces the chance of excoriation. Nonsterile, 2-ply. Individually packaged in a box of 50.

NC20421 (50)

ADHESIVE REMOVER WIPES

Clean surfaces for redressing.

Gently removes adhesive residue. Contains skin-softening aloe. Nonsterile, 2-ply. Individually packaged in a box of 100.

NC20418 (100)

KINESIO® TEX GOLD™ TAPE

High-quality, elastic tape designed for maximum range of motion.

This durable, latex-free tape is water-resistant and can be worn up to four days without irritation, even during strenuous exercise or hydrotherapy. The heat-activated adhesive will not leave any residue when removed. Kinesio® Tape has 140% elasticity, the same basic elasticity found in skin and muscles. This allows the tape to work with the body, providing full support and stability without adversely affecting the natural range of motion. Ideal for patellofemoral treatment, shin splints, plantar fasciitis, rotator cuff injuries and carpal tunnel syndrome.

Kinesio® Tex Gold™ 16.4' (5m) rolls

NC89035-202	Beige	2" (5.1cm)	(6)
NC89035-203	Blue	2" (5.1cm)	(6)
NC89035-204	Red	2" (5.1cm)	(6)
NC89035-205	Black	2" (5.1cm)	(6)

Kinesio® Tex Gold™ 103' (31.4m) Bulk roll

NC89035-332	Beige	2" (5.1cm)	(1)
-------------	-------	------------	-----



CLINICAL KINESIO® TAPING DVD

Introduction, explanation and application of the Kinesio® Taping Method.


Certified Kinesio® Taping Instructor Jim Wallis, MS, ATC, demonstrates step-by-step application of Kinesio® Taping Techniques. Color, 41 minutes. Copyright 2000. NC89048



CLINICAL THERAPEUTIC APPLICATIONS OF THE KINESIO® TAPING METHOD

Corrective application techniques for advanced Kinesio® taping.


Describes six corrective application techniques in detail: mechanical, fascia, space, ligament/tendon, functional and lymphatic. Addresses specific treatment of clinical conditions. Soft cover, 198 pages with numerous photographs. By Kenzo Kase, DC, Jim Wallis, MS, ATC and Tsuyoshi Kase. Copyright 2003. NC89047



NEW! KINESIO® TAPING IN PEDIATRICS

Fundamentals and whole body taping for infants through adolescents.


By Kenzo Kase, DC, Founder of the Kinesio® Taping Method, Patricia Martin, PT, CKTI and Audrey Yasukawa, MOT, OTR/L, CKTI. Copyright 2006. NC60012



NEW! KINESIO® TAPING FOR LYMPHOEDEMA AND CHRONIC SWELLING

Step-by-step instructions for over 108 applications.

By Kenzo Kase, DC, founder of the Kinesio® Taping Method, Kim Rock Stockheimer, MS, OTR/L, CKTI. Acknowledgement by Professor Neil Pillar. Copyright 2006. See page 85 for a more detailed description. NC60011



KINESIO® PRO SCISSORS

Great for cutting sticky items.

High-carbon, stainless steel cutting surfaces on these blades remain sharp for years. Coated with black carbon and fluorine resin that enables easy, clean cuts without leaving adhesive residue. NC89046



MEFIX®

Use as tape or to secure dressings.

Non-woven, 100% polyester fabric with a non-irritating adhesive. Latex free. Porous and non-absorbing. Minimizes fluid strike-through. Use to affix gauze, tubes, cannulas, topical dressings and drainage devices. Each roll measures 11 yds. (10m).

NC65210	1" (2.5cm)
NC65220	2" (5.1cm)
NC65240	4" (10cm)
NC65260	6" (15cm)
NC65280	8" (20cm)



FOAM UNDERWRAP

Protects skin from tape irritation.

Strong, waterproof and washable underwrap conforms to the body and accommodates tissue expansion. Non-adhesive and latex free. Each roll measures 2 3/4" x 30" (7.0cm x 76cm).

NC20700-48 Box of 48



LIGHTPLAST® PRO

Elastic tape is ideal for strapping ankles, wrists and hands.

Offers comfortable, light and cool protection. Cotton/spandex blend conforms to body movement. Resists moisture, holds securely. Each roll measures 5 yds. (4.6m).

NC12622-10	1" (2.5cm)	48 rolls
NC12622-15	1 1/2" (3.8cm)	32 rolls
NC12622-20	2" (5.1cm)	24 rolls
NC12622-30	3" (7.6cm)	16 rolls



Provides compression and support as a soft restriction splint.

No scissors required.

Make quick and easy light compression wraps.

Secure splints while providing edema compression.

Latex-Free Color Pack

DEMA™ WRAP

Lightweight, cohesive wrap provides support and compression.

Dema™ Wrap is comfortable to wear and easy to apply. Conforms to even the most difficult-to-bandage areas. Dressings are held in place without tape, pins or clips because Dema™ Wrap adheres only to itself and will not stick to hair or skin. Use to support sprains/strains and also apply compression for edema restriction. Economical for use as a bandage and as a replacement for tape when securing dressings and catheters. Rolls measure 5 yds. (4.6m) long. Also available in Latex-Free beige or Color Pack.

Standard Beige	Width	
NC68500	1" (2.5 cm)	30 Rolls
NC68501	2" (5.1 cm)	36 Rolls
NC68502	3" (7.6 cm)	24 Rolls
NC68503	4" (10 cm)	18 Rolls
Latex-Free Beige	Width	
NC68510	1" (2.5 cm)	30 Rolls
NC68511	2" (5.1 cm)	36 Rolls
Latex-Free Color Pack	Width	
NC68513	2" (5.1 cm)	36 Rolls



MOLINPIC TAPE

Tape provides excellent soft-tissue support.

Ideal for use as flexible support when weaning off a rigid splint, recovering from sports injuries or when support for weak tendons is required. Self-adherent, elastic tape won't stick to skin, clothes or hair. Available in red or blue. Each roll measures 2" x 4 yds. (5.1 cm x 3.7 m).

4 yd. (4.1 m) Rolls

NC12648-2	Red	2" (5.1 cm)
NC12649-2	Blue	2" (5.1 cm)



ECONOMY COTTON ELASTIC BANDAGES

Economical, latex-free bandages are a great value.

These latex-free, woven cotton bandages can be used for support and compression, or to secure dressings, ice, splints or casts. Wash bandage with a mild detergent and warm water on a gentle cycle. Each individually sealed bandage includes two elastic clips. Stretched length measures 4½ yds. (4.1 m).

Box of 10 Rolls - 4½ yd. (4.1 m)

NC20082	2" (5.1 cm)
NC20083	3" (7.6 cm)
NC20084	4" (10 cm)
NC20086	6" (15 cm)



PROFESSIONAL COTTON ELASTIC BANDAGES

Latex-free bandages provide controlled pressure.

Similar to the Economy Cotton Elastic Bandages sold left, these Professional bandages have about 30% less stretch. Woven cotton construction ensures improved comfort, reusability and breathability. Each individually sealed bandage includes two elastic clips. Stretched length measures 5 yds. (4.6 m).

Box of 10 Rolls - 5 yd. (4.6 m)

NC20072	2" (5.1 cm)
NC20073	3" (7.6 cm)
NC20074	4" (10 cm)
NC20076	6" (15 cm)



FABRIFOAM® NuStimWrap™

Four-way stretch wrap offers support and compression.

Made of a patented fabric and foam composite material that provides moderate elasticity. Use these non-neoprene, latex-free wraps to hold hot and cold packs. Sold in boxes of three wraps.

Box of 3 Wraps

NC20092-1	2½" x 18" (6.4 x 46 cm)
NC20092-2	2½" x 36" (6.4 x 91 cm)
NC20092-3	2½" x 48" (6.4 x 120 cm)
NC20092-4	4" x 18" (10 x 46 cm)
NC20092-5	4" x 36" (10 x 91 cm)
NC20092-6	4" x 48" (10 x 120 cm)
NC20092-7	6" x 18" (15 x 46 cm)
NC20092-8	6" x 36" (15 x 91 cm)
NC20092-9	6" x 48" (15 x 120 cm)




NEW! IMAK® ARTHRITIS GLOVES

Mild compression helps control and decrease joint pain and swelling.

Use IMAK® Arthritis Gloves to relieve the aches, pains and stiffness of arthritis. May also help with hand fatigue and other types of hand and wrist pain, neuropathy of the hands and poor circulation. Made from soft, breathable cotton LYCRA®, which allows skin to breathe and keeps hands dry. Open fingertips permit better gripping and tactile sensation. Comfortable enough to use all day and night. Sold in pairs. Hand wash and air dry. Latex free.

		Width at Knuckles
NC53108-1	Small	Up to 3 1/8" (Up to 7.9cm)
NC53108-2	Medium	Up to 3 1/2" (Up to 8.9cm)
NC53108-3	Large	Up to 4" (Up to 10cm)



IMAK® Arthritis Glove Sizing

To size, measure hand width at knuckles as indicated.



DEMA® GRIP COMPRESSION STOCKINETTE

Economical cotton compression stockinette for edema care.

Made of cotton and elastic stockinette. Latex free. Dema® Grip self-adjusts to provide consistent pressure over contours. Can be washed or autoclaved. Each roll measures 33' (10m). Available in natural and beige.

Natural	Size	Flat Width	Beige	Size	Flat Width
NC42500	A	2" (5.1cm)	NC42601	B	2 1/2" (6.4cm)
NC42501	B	2 1/2" (6.4cm)	NC42602	C	2 3/4" (6.9cm)
NC42502	C	2 3/4" (6.9cm)	NC42603	D	3" (7.6cm)
NC42503	D	3" (7.6cm)	NC42604	E	3 1/2" (8.9cm)
NC42504	E	3 1/2" (8.9cm)	NC42605	F	4" (10cm)
NC42505	F	4" (10cm)	NC42607	G	4 3/4" (12cm)
NC42507	G	4 3/4" (12cm)			
NC42508	J	6 3/4" (17cm)			
NC42509	K	8" (20cm)			
NC42510	L	10" (25cm)			
NC42511	M	12 3/4" (30cm)			

Using the flat width to determine size:


- Measure circumference of the area needing compression.
- Divide measurement in half.
- Select the size closest to the measurement.

IMPROVED! ISOTONER® THERAPEUTIC GLOVES

Smooth, two-way stretch fabric provides gentle, even compression to help prevent edema.

The Arthritis Foundation® found that "the gloves potentially can provide an improvement in ease of use as an addition to other treatments in edema control and grip strength." Made of 80% nylon and 20% spandex. Latex free. Available in either **Open** or **Full Finger** styles. Open finger design offers extra mobility. All seams are on the outside of the glove for smooth comfort. Area around base of the thumb is free of seams. Provides moderate compression of approximately 25 to 35mm Hg pressure. Compression is estimated from palm circumference at MP joints. Machine wash in cold water on delicate cycle. Lay flat to dry. Gloves are unisex sized for both men and women. Sold in pairs. Guaranteed for one year against defects in material and workmanship.

Open Finger		MP Circumference
NC53022-0	X-Small	4 1/2" to 5 1/2" (11 to 14cm)
NC53022-1	Small	5 1/2" to 6 1/2" (14 to 17cm)
NC53022-2	Medium	6 1/2" to 7 1/2" (17 to 19cm)
NC53022-3	Large	7 1/2" to 8 1/2" (19 to 22cm)
Full Finger		MP Circumference
NC53023-0	X-Small	4 1/2" to 5 1/2" (11 to 14cm)
NC53023-1	Small	5 1/2" to 6 1/2" (14 to 17cm)
NC53023-2	Medium	6 1/2" to 7 1/2" (17 to 19cm)
NC53023-3	Large	7 1/2" to 8 1/2" (19 to 22cm)



Isotoner® Glove Sizing

Measure circumference around the MP joints as indicated to determine size.



Isotoner® Full Finger Gloves



Isotoner® Open Finger Gloves

Features all outside seams and finished edges. Base of thumb is seam free.

Open cuff for easy donning and doffing.

Open fingers help improve mobility.

Material extends to cover painful joints.

Comfortable glove fits
left or right hand.



NEW! ELVAREX® SEAMLESS GLOVE

Seamless glove delivers increased comfort.

Seamless design provides more comfort in the management of lymphedema. Features include finished finger stub ends for increased durability, correct anatomical design, reinforced zone between the thumb and index finger, and reinforced wristband for easier donning and doffing. Reduced compression at the wrist prevents a doubling of compression when worn with a JOBST® Ready-to-Wear Armsleeve (page 83). Compression also is reduced at the finger stubs, preventing a tourniquet effect. Glove has a mean compression average of 20-30 mmHg. Sold individually. To size, measure palm and wrist circumference. Beige color.

	Size	Palm Circumference	Wrist Circumference
NC60061-1	1	6¼" to 6⅞" (16 to 17 cm)	5" to 5¾" (13 to 15 cm)
NC60061-2	2	6⅞" to 7½" (17 to 19 cm)	5¾" to 6¼" (15 to 16 cm)
NC60061-3	3	7½" to 8" (19 to 20 cm)	6¼" to 6⅞" (16 to 17 cm)
NC60061-4	4	8" to 8⅞" (20 to 22 cm)	6⅞" to 7½" (17 to 19 cm)
NC60061-5	5	8⅞" to 9¼" (22 to 23 cm)	7½" to 8" (19 to 20 cm)
NC60061-6	6	9¼" to 9¾" (23 to 25 cm)	8" to 8⅞" (20 to 22 cm)

SWELL SPOTS™

These inserts are ideal for treating swollen and fibrotic areas that are not helped with simple bandaging.

Swell Spots™ feature an enhanced, quilted channel system. Made of cotton-based fabric and soft, lightweight foam chips. Use to treat lymphedema and tissue fibrosis. Save time, no need to make chip bags. Machine washable and latex free. Patent Number 5,976,099.

Dorsal Spot: Use primarily on the dorsum of the hand and foot. Works well under edema gloves and compression wraps. Small Dorsal Spot measures 3½" x 3⅞" (10 x 7.9 cm). Large Dorsal Spot measures 5⅞" x 3⅞" (15 x 9.0 cm).

Oval Spot: Ideal for areas such as the axilla and inside a carpal tunnel splint. Measures 7½" x 3⅞" (19 x 7.9 cm).

Katina Spot: Ideal for areas such as the ulna or tibia. Small Katina Spot measures 15" x 3⅞" (38 x 10 cm). Large Katina Spot measures 16½" x 7½" (42 x 18 cm).

- NC60020-2 Small Dorsal Spot
- NC60020-3 Large Dorsal Spot
- NC60020-6 Oval Spot
- NC60020-8 Small Katina Spot
- NC60020-7 Large Katina Spot



COMPRESSION GLOVES

Made of unbleached cotton, nylon and spandex. Provides gentle compression (approximately 15-25 mmHg pressure) to relieve lymphedema and other types of hand edema. Frequent machine washing and drying will enhance the performance of these gloves. All glove styles are open tipped. Sold individually. To size, measure both MP and wrist circumferences. Refer to the sizing charts below.

Compression Gloves Sizing Chart

Sizing for Style 02			Sizing for Styles 14 and 34		
Size	Fits MPs	Fits Wrists	Size	Fits MPs	Fits Wrists
S	17.5 to 18.4 cm (6⅞" to 7¼")	14.0 to 15.0 cm (5½" to 5⅞")	S	18.4 to 19.5 cm (7¼" to 7⅞")	13.5 to 14.5 cm (5⅜" to 5⅞")
M	18.4 to 19.5 cm (7¼" to 7⅞")	14.5 to 15.6 cm (5⅞" to 6¼")	M	19.5 to 20.5 cm (7⅞" to 8⅞")	14.0 to 14.9 cm (5½" to 5⅞")
L	19.5 to 20.5 cm (7⅞" to 8⅞")	15.0 to 16.0 cm (5⅞" to 6¼")	L	20.5 to 21.4 cm (8⅞" to 8⅞")	14.4 to 15.6 cm (5⅞" to 6¼")
XL	20.5 to 21.5 cm (8⅞" to 8⅞")	15.6 to 16.5 cm (6⅞" to 6½")	XL	21.4 to 22.5 cm (8⅞" to 8⅞")	15.0 to 16.0 cm (5⅞" to 6¼")
XXL	21.5 to 22.5 cm (8⅞" to 8⅞")	16.0 to 17.0 cm (6⅞" to 6⅞")	XXL	22.5 to 23.5 cm (8⅞" to 9¼")	15.6 to 16.5 cm (6⅞" to 6½")

Regular Length:*

Style 02: Outside seams. Long finger length to nail area. Overall length to wrist crease. Darts on dorsum for extra compression.

*Also available with the finger length to PIP. Call for information.



Style 14: Outside seams. Long finger length to nail area beyond DIP. Overall length to wrist crease. No darts on dorsum side. Fits left or right hand.

*Also available with the finger length to PIP. Call for information.

Regular Length Style 02

Left	Right	
NC60030-1	NC60030-2	S
NC60030-3	NC60030-4	M
NC60030-5	NC60030-6	L
NC60030-7	NC60030-8	XL
NC60030-9	NC60030-10	XXL

Regular Length Style 14 Fits left or right hand.

NC60031-1	S
NC60031-2	M
NC60031-3	L
NC60031-4	XL
NC60031-5	XXL

Long Length:*

Style 34: Glove length is 6.0 to 7.0 cm (2⅞" to 2¾") proximal to wrist. Outside seams, long finger length to nail area. No darts on dorsum side. Fits left or right hand.

*Also available with the finger length to PIP. Darts on dorsum for extra compression. Fits left or right hand. Call for information.

Long Length Style 34

Fits left or right hand.

NC60032-1	Small
NC60032-2	Medium
NC60032-3	Large
NC60032-4	X-Large
NC60032-5	XX-Large



Components	Leg Kit	Arm Kit	Lower Leg Kit
Tricofix® D5		(1) 6cm x 20m (2.4" x 22 yd)	
Tricofix® E6			(1) 8cm x 20m (3.2" x 22 yd)
Tricofix® F7	(1) 10cm x 20m (3.9" x 22 yd)		
Elastomull®	(20) 4cm x 3.7m (1.6" x 4 yd)	(20) 4cm x 3.7m (1.6" x 4 yd)	(20) 4cm x 3.7m (1.6" x 4 yd)
Artiflex®	(2) 10cm x 3.4m (4.0" x 3.7 yd) (2) 15cm x 3.4m (6.0" x 3.7 yd)	(2) 10cm x 3.4m (4.0" x 3.7 yd) (2) 15cm x 3.4m (6.0" x 3.7 yd)	(3) 10cm x 3.4m (4.0" x 3.7 yd)
Comprilan®	(1) 6cm x 5m (2.4" x 5.5 yd) (1) 8cm x 5m (3.2" x 5.5 yd) (3) 10cm x 5m (4.0" x 5.5 yd) (4) 12cm x 5m (4.7" x 5.5 yd)	(1) 6cm x 5m (2.4" x 5.5 yd) (1) 8cm x 5m (3.2" x 5.5 yd) (5) 10cm x 5m (4.0" x 5.5 yd)	(1) 6cm x 5m (2.4" x 5.5 yd) (1) 8cm x 5m (3.2" x 5.5 yd) (3) 10cm x 5m (4.0" x 5.5 yd)
Leukoplast® hospital Surgical Tape	(1) 1.25cm x 9.2m (0.5" x 10 yd)	(1) 1.25cm x 9.2m (0.5" x 10 yd)	(1) 1.25cm x 9.2m (0.5" x 10 yd)
JOBST® Foam Rubber Pad	(1) Kidney Shape, 9cm x 1cm (3.5" x 0.4")	(1) Kidney Shape, 9cm x 1cm (3.5" x 0.4")	(1) Kidney Shape, 9cm x 1cm (3.5" x 0.4")
	NC60001	NC60002	NC60003



NEW! LymphCare Kits

Kits offer a complete range of bandage products for lymphedema management.

LymphCare Kits offer a complete solution for compression bandaging of upper and lower extremities. Three kits are available: the **Leg Kit**, the **Arm Kit**, and the **Lower Leg Kit**. Components are listed in chart on left.

Components	Set 4 - Large Arms	Set 5 - Small Legs	Set 6 - Large Legs
tg® Tubular Bandage	(1) Size 7 (7.0cm x 20m)	(1) Size 9 (8.6cm x 20m)	(1) Size 9 (8.6cm x 20m)
Mollelast® Bandage	(20) 4cm x 4m (1.6" x 4.4 yd)	(20) 4cm x 4m (1.6" x 4.4 yd)	(20) 6cm x 4m (2.4" x 4.4 yd)
Kompres® Foam Rubber Bandage	(1) 10cm x 1m (4.0" x 1.1 yd) 1 cm thickness (0.4")	(1) 8cm x 1m (3.2" x 1.1 yd) 0.5 cm thickness (0.2")	(1) 10cm x 1m (4.0" x 1.1 yd) 1 cm thickness (0.4")
Cellona® Synthetic Padding	(4) 6cm x 3m (2.4" x 3.3 yd) (4) 10cm x 3m (4.0" x 3.3 yd)		
Rosidal® soft		(2) 10cm x 2.5m (4.0" x 2.7 yd) 0.3 cm thickness (.11") (2) 15cm x 2.5m (6.0" x 2.7 yd) 0.4 cm thickness (.16")	(3) 10cm x 2.5m (4.0" x 2.7 yd) 0.3 cm thickness (.11") (3) 15cm x 2.5m (6.0" x 2.7 yd) 0.4 cm thickness (.16")
Rosidal® K Bandage	(2) 6cm x 5m (2.4" x 5.5 yd) (2) 10cm x 5m (4.0" x 5.5 yd) (3) 12cm x 5m (4.7" x 5.5 yd)	(2) 6cm x 5m (2.4" x 5.5 yd) (2) 8cm x 5m (3.2" x 5.5 yd) (2) 10cm x 5m (4.0" x 5.5 yd) (4) 12cm x 5m (4.7" x 5.5 yd)	(2) 6cm x 5m (2.4" x 5.5 yd) (2) 8cm x 5m (3.2" x 5.5 yd) (4) 10cm x 5m (4.0" x 5.5 yd) (6) 12cm x 5m (4.7" x 5.5 yd)
Idealbinde®		(1) 20cm x 5m (7.9" x 5.5 yd)	(2) 20cm x 5m (7.9" x 5.5 yd)
	NC60005	NC60006	NC60007



NEW! ROSIDAL® LYPHSET

Ready-to-use, time-saving sets.

These comprehensive sets contain all of the necessary components for a professional lymphological compression dressing. Three sets are available: Set 4 (for large arms), Set 5 (for small legs), and Set 6 (for large legs).

NEW! JOBST® READY-TO-WEAR ARMSLEEVE

Exceptionally soft and comfortable.

This durable, lightweight knit armsleeve assists in managing hypertrophic scarring and conditions that require compression following trauma, surgery or radiation therapy. Available in two compressions. The 15-20 mmHG line is for patients who cannot tolerate higher compression. Hand washable. Beige. To size, measure circumference at wrist and arm, as indicated below.



	Compression	Wrist	Mid-Lower	Mid-Upper
NC60064-1 S	15-20 mmHg	5 3/8" to 7 3/8"	7" to 9 1/8"	8 1/2" to 15 1/2"
NC60064-2 M	15-20 mmHg	6 1/8" to 8 1/2"	9 1/4" to 11 1/4"	10 3/8" to 17 3/4"
NC60064-3 L	15-20 mmHg	6 7/8" to 9 1/4"	11 3/8" to 13 1/2"	12 3/4" to 19 7/8"
	Compression	Wrist	Mid-Lower	Mid-Upper
NC60065-1 S	20-30 mmHg	5 3/8" to 7 3/8"	7" to 9 1/8"	8 1/2" to 15 1/2"
NC60065-2 M	20-30 mmHg	6 1/8" to 8 1/2"	9 1/4" to 11 1/4"	10 3/8" to 17 3/4"
NC60065-3 L	20-30 mmHg	6 7/8" to 9 1/4"	11 3/8" to 13 1/2"	12 3/4" to 19 7/8"

NEW! JOBST® READY-TO-WEAR GAUNTLET

Reduced compression at the wrist.

This soft, lightweight knit gauntlet is appropriate for edema and mild to moderate lymphedema. Provides a seamless fit when worn with the Ready-to-Wear Armsleeve (see left). May be worn alone. Available in 15-20 mmHG and 20-30 mmHG. Hand washable. Beige. To size, measure circumference at palm and wrist, as indicated below.



	Compression	Palm	Wrist
NC60062-1 S	15-20 mmHg	5 7/8" to 7 3/8"	5 3/8" to 7 3/8"
NC60062-2 M	15-20 mmHg	7 1/2" to 8 5/8"	6 1/8" to 8 1/2"
NC60062-3 L	15-20 mmHg	8 3/4" to 10 1/4"	6 7/8" to 9 1/4"
	Compression	Palm	Wrist
NC60063-1 S	20-30 mmHg	5 7/8" to 7 3/8"	5 3/8" to 7 3/8"
NC60063-2 M	20-30 mmHg	7 1/2" to 8 5/8"	6 1/8" to 8 1/2"
NC60063-3 L	20-30 mmHg	8 3/4" to 10 1/4"	6 7/8" to 9 1/4"

Reduced compression at the wrist for a seamless fit with armsleeves.

**TG-TUBULAR™
BANDAGE
STOCKINETTE**

A great protective cover under all casting/wrapping applications, from small fingers to large body casts.

Fine, seamless, knitted tube made of 67% bleached cotton and 33% viscose. Material is tightly meshed and can be stretched laterally without creases.

For fingers, toes and children's hands:
Size 5: 5.6 cm x 20 m (2 3/16" x 22 yds.)

NC60024-5 (1)
NC60024-5C (10)

For arms, feet, and children's legs:
Size 6: 6.4 cm x 20 m (2 1/2" x 22 yds.)

NC60024-6 (1)
NC60024-6C (10)

For legs, thighs and arms:
Size 7: 7.0 cm x 20 m (2 3/4" x 22 yds.)

NC60024-7 (1)
NC60024-7C (18)

For thighs and heads:
Size 9: 8.5 cm x 20 m (3 3/8" x 22 yds.)
NC60024-9 (1)

**KOMPRES® FOAM RUBBER**

Porous, synthetic, latex foam rubber is permeable to air and water. Lightweight with a soft surface. Cut to create customized shapes. Autoclavable at 250°F (121°C).

Roll:
8.0 cm x 2.0 m x 5.0 mm (3 1/8" x 6 1/2" x 3/16")

NC60029-1 (1)
NC60029-1C (30)

Kidney shape: Size 0
10 mm x 9.0 cm x 5.0 cm (7/16" x 3 1/2" x 2")

NC60029-2 (1)
NC60029-2C (75)

Kidney shape: Size 1
10 mm x 12 cm x 6.5 cm (7/16" x 12" x 6 1/2")

NC60029-3 (1)
NC60029-3C (30)

Sheet:
10 mm x 1.0 m x 50 cm (7/16" x 1.1 yd. x 20")

NC60029-4 (1)
NC60029-4C (10)

**ELASTOMULL®**

This flexible, comfortable, elastic bandage stays in place even with constant movement. 100% cotton, seamless, knitted gauze. High-density fabric has a soft, rich texture. Stretches to approximately twice its size. Non-sterile.

2.5 cm x 3.7 m (1" x 12')

NC20290-1 (24)
NC20290-1C (192)

5.1 cm x 3.7 m (2" x 12')

NC20290-2 (12)
NC20290-2C (96)

7.6 cm x 3.7 m (3" x 12')

NC20290-3 (12)
NC20290-3C (96)

10 cm x 3.7 m (4" x 12')

NC20290-4 (12)
NC20290-4C (96)

**TRICOFIX®**

This 100% cotton, seamless, knitted tube gauze will not fray when cut. Stretches lengthwise and crosswise for easy application. Sold in 20 m (22 yd) lengths.

6.0 cm x 20 m (2 3/8" x 22 yds.)

NC20259-6 (1)
NC20259-6C (10)

8.0 cm x 20 m (3 1/8" x 22 yds.)

NC20259-8 (1)
NC20259-8C (10)

10 cm x 20 m (4" x 22 yds.)

NC20259-10 (1)
NC20259-10C (10)

12 cm x 20 m (4 3/4" x 22 yds.)

NC20259-12 (1)
NC20259-12C (10)

**ROSIDAL™ SOFT**

Ideally suited for padding under lymphological bandages. The interlocking open pores help ensure even pressure distribution of short-stretch bandages. Polyurethane foam is highly breathable and permeable to moisture. Tear resistant. Does not lose its elasticity after repeated washing. Latex free.

10 cm x 3 mm x 2.5 m (3 15/16" x 1/8" x 2 3/4 yds.)

NC60040-1 (1)
NC60040-1C (30)

10 cm x 4 mm x 2.5 m (3 15/16" x 3/16" x 2 3/4 yds.)

NC60040-2 (1)
NC60040-2C (24)

12 cm x 4 mm x 2.5 m (4 3/4" x 3/16" x 2 3/4 yds.)

NC60040-3 (1)
NC60040-3C (16)

15 cm x 4 mm x 2.5 m (5 7/8" x 3/16" x 2 3/4 yds.)

NC60040-4 (1)
NC60040-4C (14)

**CELLONA®
SYNTHETIC PADDING**

Use this self-conforming padding under any bandaging medium. Water-repellent, non-woven, polyester fabric. Each roll measures 3.7 m (4 yds.) long.

5.1 cm x 3.7 m (2" x 4 yds.)

NC60027-1 (12)
NC60027-1C (72)

7.6 cm x 3.7 m (3" x 4 yds.)

NC60027-2 (12)
NC60027-2C (72)

**ARTIFLEX®**

Use as padding under bandages to ensure even distribution of pressure. Non-absorbent. Autoclave at 250°F (121°C) for sterile application. Made of non-woven 40% polyester, 40% polypropylene and 20% polyethylene. Roll measures 10 cm x 3.4 m (3 15/16" x 3 3/4 yds.).

NC60028-1 (1)
NC60028-30 (30)

**NEW! JOBST® FOAM PADS
AND BANDAGES**

These synthetic, latex foam pads and bandages distribute pressure evenly. Reusable. Hand wash and air dry.

Foam Pad: Kidney shape
9.0 x 1.0 cm (3 1/2" x 2/5")

NC60068 (1)

Foam Bandage Roll

10 cm x 1.0 m x 1.0 cm (3 15/16" x 1.1 yd. x 3/8")
NC60069 (1)





ROSIDAL® K BANDAGE

This durable bandage wrap is made of 100% cotton. One-way, medium stretch, longitudinally elastic with 90% maximum stretch. Two clips included. Machine wash. Each roll measures 5.0m (16½') long.

6.0 cm (2⅜")

NC60022-1 (1)

NC60022-1C (108)

8.0 cm (3⅞")

NC60022-2 (1)

NC60022-2C (108)

10 cm (3⅞")

NC60022-3 (1)

NC60022-3C (90)

12 cm (4¾")

NC60022-4 (1)

NC60022-4C (54)



ISOBAND® FINISHED EDGE BANDAGE

Large wrap for use with chip bags. One-way short stretch, longitudinally elastic with 90% maximum stretch for lasting comfort. Porous weave allows skin to breathe. Made of 85% cotton and 15% cellulose for superior comfort. May be washed and dried flat without losing any of its properties. Each roll measures 5.0m (16½') long.

14 cm (5½")

NC60023-1 (1)

NC60023-1C (50)

20 cm (7⅞")

NC60023-2 (1)

NC60023-2C (35)

NEW! KINESIO® TAPING FOR LYMPHOEDEMA AND CHRONIC SWELLING

Introduces unique techniques that compliment other "standard" lymphoedema treatments.


This easy-to-read book demonstrates step-by-step instructions of over 108 common and not-so-common applications associated with lymphoedema, chronic swelling and scar management. Kinesio® Tex Tape is applied directly to the skin and assists in the movement of muscles, repositioning of joints and alleviation of pain.

The Kinesio® Taping Method assists in the "lifting or creation of space" within the superficial skin, which decreases pressure while opening initial lymphatics. By allowing maximum muscle contraction and relaxation, efficiency of the deeper lymphatics also may be improved.

Full-color, spiral bound, 172 pages. Written by Kenzo Kase, DC and Kim Rock Stockheimer, MS, OTR/L, CKTI. Acknowledgement by Professor Neil Piller. Copyright 2006.

NC60011

See page 79 for more Kinesio® products.



KOMPRIBAND™

Prolonged memory offers excellent conformity. One-way, short stretch, medium compression wrap. Made of 100% unbleached cotton. Machine wash and dry. Each roll measures 5.0m (16½') long.

NC60016-1 6.0 cm (2⅜")

NC60016-2 8.0 cm (3⅞")

NC60016-3 10 cm (3⅞")

NC60016-4 12 cm (4¾")



FINGERBAND™, TOEBAND™ AND HANDBAND™

Ideal for wrapping infant extremities. One-way, short stretch, medium compression wrap. Made of 100% unbleached cotton. Machine wash and dry. Rolls are 3.5m (11½') long.

FingerBAND™

NC60017-1 3.0 cm (1⅜")

ToeBAND™

NC60017-2 19 mm (¾")

HandBAND™

NC60017-3 4.0 cm (1⅞")



IDEALBINDE®

One-way, short stretch, strong compression wrap with 90% extensibility, durable elasticity. Ideal for use with chip bags. Made of 86% cotton and 14% viscose. Machine wash. May be boiled or sterilized at 250° F (121°C). Each roll measures 15 cm x 5.0m (5⅞" x 16½') long.

NC60021 (1)

NC60021C (60)



COMPRILAN®

Open weave, one-way, short stretch wrap offers consistent compression while allowing breathability. Longitudinally elastic with 90% maximum stretch. Durable, easy-care, 100% cotton. Machine wash and dry. Each roll measures 5.0m (16½') long.

4.0 cm (1½")

NC20265 (1)

NC20265C (20)

6.0 cm (2⅜")

NC20261 (1)

NC20261C (20)

8.0 cm (3⅞")

NC20263 (1)

NC20263C (20)

10 cm (3⅞")

NC20260 (1)

NC20260C (20)

12 cm (4¾")

NC20264 (1)

NC20264C (20)



NEW! X-STATIC® TUBULAR COMPONENT

Forms a non-compressive protective layer between skin and bandages or casts.


Made of 10% X-Static® and 90% cotton or polyester. X-Static® fibers are 99.9% pure elemental silver which offers natural antimicrobial protection to help eliminate bacteria and odors. Available in 25 yd. (23m) rolls, in four widths. Latex free.

NC60046-2 2" (5.1 cm)

NC60046-4 4" (10 cm)

NC60046-3 3" (7.6 cm)

NC60046-6 6" (15 cm)





TENDERWRAP® UNNA BOOT BANDAGE

Non-hardening, zinc oxide formulation.

Soft, non-raveling, 100% cotton gauze is impregnated with a non-hardening zinc oxide paste. Bandage remains flexible, encouraging patient ambulation and compliance. Conforms well to body contours and provides even compression. Gelatin-free and preservative-free bandages are individually wrapped in non-sterile packaging. Sold in boxes of 12.

3" x 10 yds. (7.6cm x 9.1m)

NC95209-3 (12)

4" x 10 yds. (10cm x 9.1m)

NC95209-4 (12)



NEW! LISTER BANDAGE SCISSORS

Two sizes to suit your needs.

Blunt-sided utility scissors are made of stain-resistant steel to resist corrosion and rusting. Ideal for removing gauze or bandages. Blunt end decreases the chance of accidental skin injuries. Autoclavable. Total length of the scissors measures 5½" or 7¼" (14 or 18cm). For right hand use.

NC12742 5½" (14cm)

NC12742-1 7¼" (18cm)

HIBICLENS® ANTISEPTIC/ANTIMICROBIAL SKIN CLEANSER

Reduce the risk of cross infection.

Hibiclens® can be used as a surgical hand scrub, a healthcare personnel handwash, a patient pre-operative skin preparation, and as a skin wound cleanser and general skin cleaner. Hibiclens® is an effective defense against the spread of MRSA and other staph infections. For external use only. Do not use near eyes, ears or mouth.

NC12690 4 fl. oz. (118ml)

NC12690-1 32 fl. oz. (946ml)



CURI-STRIP® ADHESIVE WOUND CLOSURES

Hypoallergenic adhesive secures wounds.

Conformable, spunbonded nylon backing that is free of reinforcing filaments for better movement. Helps minimize premature "pop-off" and formation of stress blisters. Highly porous backing for increased oxygen exposure. Superior moisture vapor transmission for minimal skin maceration.

5 per envelope, 50 envelopes per box

NC95200-1 ⅛" x 3" (3.2mm x 7.6cm) (250)

6 per envelope, 50 envelopes per box

NC95200-2 ¼" x 1½" (6.4mm x 3.8cm) (300)

NC95200-5 ½" x 4" (13mm x 10cm) (300)

Bulk box

NC95200-6 1" x 4" (2.5 x 10cm) (100)



NEW! SUTURE REMOVAL KIT

Disposable kits ensure sterility.

Sterile kit contains metal-insert forceps, Littauer Scissors and a cotton-filled gauze sponge. Ideal for wound care following injury or surgery. Designed for one-time use.

NC20149 Single kit

NC20150 Box of 50 kits



EXAM GLOVES

Available in a variety of choices.

Sterile Gloves are individually wrapped and sold in pairs.

Sterile: Box of 50

NC20511-0 Small

NC20511-1 Medium

NC20511-2 Large

Sterile, Latex-Free: Box of 50

NC20572-1 Small

NC20572-2 Medium

NC20572-3 Large

Nitrile Non-Sterile, Latex-Free and Powder-Free: Box of 100

NC20510-1 Small

NC20510-2 Medium

NC20510-3 Large

NC20510-4 X-Large

Non-Sterile, Latex-Free: Box of 100

NC20517-1 Small

NC20517-2 Medium

NC20517-3 Large

NC20517-4 X-Large

Non-Sterile, Latex: Box of 100

NC20518-1 Small

NC20518-2 Medium

NC20518-3 Large

NC20518-4 X-Large

Non-Sterile and Powder-Free: Box of 100

NC20519-1 Small

NC20519-2 Medium

NC20519-3 Large

NC20519-4 X-Large

Non-Sterile Vinyl, Powder-Free: Box of 100

NC20520-1 Small

NC20520-2 Medium

NC20520-3 Large

NC20520-4 X-Large



SKIN PREP/BARRIER WIPES AND ADHESIVE REMOVER WIPES

Protect and prepare skin for dressing.

Skin Prep Wipes protect the skin; Adhesive Remover Wipes remove residue from bandages and tapes. For more information, see page 78 or visit www.ncmedical.com.



KERLIX® ROLLS

Unique soft gauze cushions and conforms.

Ideal for bandaging heads, limbs and difficult-to-dress wounds such as burns, plastic and orthopedic wounds. Made of prewashed, fluff-dried, 100% woven gauze with unique crinkle weave pattern for loft and bulk to cushion and protect wound areas. Six-ply construction with finished edges reduces loose ends and lint. Sterile rolls are individually wrapped.

2 1/4" x 3 yds. (5.7 cm x 2.7 m)

NC20131 (12)

NC20131C (96)

4 1/2" x 4 yds. (11 cm x 3.7 m)

NC20132 (10)

NC20132C (100)



CURASORB™ CALCIUM ALGinate®

Absorbs approximately 20 times the dressing's weight.

These dressings absorb exudate from the wound and create a protective cushioning gel and optimal moist healing environment. Protective gel binds exudate and helps prevent skin maceration. Easy, painless dressing removal minimizes tissue trauma. Individually wrapped in sterile packaging.

NC95206-1 12" (30 cm) rope (5)

NC95206-2 2" x 2" (5.1 x 5.1 cm) (10)

NC95206-3 4" x 4" (10 x 10 cm) (10)



XEROFORM® PETROLATUM DRESSING

Clings and conforms to body contours.

Contains non-adherent, 3% Bismuth Tribromophenate in a non-medicinal, petrolatum blend on fine mesh gauze. Bacteriostatic action. Prevents gauze from adhering to wounds. Individually wrapped in sterile packaging.

Box of 50

NC20141 1" x 8" (2.5 x 20 cm)

NC20142 5" x 9" (13 x 23 cm)



SPONGES

Versatile sponges offer maximum absorption.

Cotton-filled sponges can be used for a wide range of dressing, cleaning and prepping procedures. Cut edges are folded, no raw edges are exposed. Non-sterile sponges are available in two styles: All-gauze and Non-woven. Variety of standard sizes and plys.

All-gauze, non-sterile sponges: Box of 200

NC20270-208 2" x 2" (5.1 x 5.1 cm) 8-ply

NC20270-308 3" x 3" (7.6 x 7.6 cm) 8-ply

NC20270-408 4" x 4" (10 x 10 cm) 8-ply

NC20270-212 2" x 2" (5.1 x 5.1 cm) 12-ply

NC20270-312 3" x 3" (7.6 x 7.6 cm) 12-ply

NC20270-412 4" x 4" (10 x 10 cm) 12-ply

Non-woven, non-sterile sponges: Box of 200

NC20273-204 2" x 2" (5.1 x 5.1 cm) 4-ply

NC20273-304 3" x 3" (7.6 x 7.6 cm) 4-ply

NC20273-404 4" x 4" (10 x 10 cm) 4-ply

All-gauze, sterile sponges: Box of 25

NC20166 3" x 3" (7.6 x 7.6 cm) 12-ply

NC20167 4" x 4" (10 x 10 cm) 12-ply

TELEFA® ISLAND DRESSINGS

Convenient, all-in-one dressing.

Ideal for use on lightly draining wounds. The soft, non-woven backing conforms around the wound and seals on all four sides to provide protection from the external environment. Non-adherent TELEFA® will not disrupt healing tissue. A secure dressing for central and peripheral IV sites. Individually wrapped in sterile packaging.

NC95202-2 4" x 5" (10 x 13 cm) (25)

NC95202-3 4" x 8" (13 x 20 cm) (25)



CURITY® NON-ADHERING DRESSINGS

Allows drainage of exudate away from wound areas.

Ideal for donor sites, delicate suture lines, graft sites and draining wounds. Oil-emulsion blend impregnated into an open-mesh, knitted fabric dressing. Individually wrapped, sterile dressings are available in two sizes.

Box of 50

NC95213-2 3" x 3" (7.6 x 7.6 cm)

Box of 24

NC95213-3 3" x 8" (7.6 x 20 cm)



CONSTANT-CLENS® DERMAL WOUND CLEANSER

Reduce the risk of cross infection.

Cleanser removes and softens necrotic tissue and debris. Isotonic and pH balanced, non-toxic and non-irritating. Contains Polysorbate 20.

8 fl. oz (237 ml) Spray Bottle

NC95201-1 (1)

NC95201-12 (12)



pH-balanced

N-TERFACE

Ideal for post-operative wound care.

Lightweight and porous, this second-skin material prevents dressing from adhering to wounds. Light tack keeps the material in position but will not pull on skin or wounds when removed. N-terface can be repositioned as needed. Individually wrapped, sterile dressings measure 4" x 4" (10 x 10 cm).

Box of 25

NC20009 4" x 4" (10 x 10 cm)

