

The four mat platforms shown here feature:

Quality Construction

- Solid oak legs (except Economy).
- Frame strengthened by steel angles.
- Standard 18" (46cm) height.
- Custom bumper strip that protects the vinyl top.



CLINICA™ ECONOMY MAT PLATFORM

An affordable choice and a great value.

- Weight capacity is 650 lbs. (295 kg).
- Call for freight quote.

| | Width x Length |
|-----------|----------------------|
| NC9200-47 | 4' x 7' (1.2 x 2.1m) |
| NC9200-57 | 5' x 7' (1.5 x 2.1m) |
| NC9200-68 | 6' x 8' (1.8 x 2.4m) |



CLINICA™ MAT PLATFORM WITH REMOVABLE MAT

Mat attaches easily with Velcro® strips.

- Weight capacity is 650 lbs. (295 kg).
- Call for freight quote.

| | Width x Length |
|-----------|----------------------|
| NC9204-47 | 4' x 7' (1.2 x 2.1m) |
| NC9204-57 | 5' x 7' (1.5 x 2.1m) |
| NC9204-68 | 6' x 8' (1.8 x 2.4m) |



CLINICA™ HYDRAULIC HAND-CRANK MAT PLATFORM

Easy-to-use hand crank folds out of the way.

- Heavy-duty, bolted legs with levelers and four corner lifting cylinders.
- Height adjusts from 23" to 35" (58 to 89cm).
- Premium, stain-resistant vinyl upholstery with knit-backing. Choice of eight colors (see page 20).
- Weight capacity is 750 lbs. (340 kg).
- Call for freight quote.



| | Width x Length |
|-----------|----------------------|
| NC9210-47 | 4' x 7' (1.2 x 2.1m) |
| NC9210-57 | 5' x 7' (1.5 x 2.1m) |
| NC9210-68 | 6' x 8' (1.8 x 2.4m) |

Premium Upholstery

- Ultra-firm, 2" (5.1cm) padding.
- Premium, scratch- and stain-resistant vinyl upholstery with knit-backing. Choice of eight colors (see page 20).



CLINICA™ UPHOLSTERED MAT PLATFORM

Six bolted legs provide excellent stability.

- Weight capacity is 650 lbs. (295 kg).
- Call for freight quote.

| | Width x Length |
|-----------|----------------------|
| NC9202-47 | 4' x 7' (1.2 x 2.1m) |
| NC9202-57 | 5' x 7' (1.5 x 2.1m) |
| NC9202-68 | 6' x 8' (1.8 x 2.4m) |



CLINICA™ FOLDING MAT PLATFORM

Ideal where space is limited.

- Projects 8" (23cm) when folded against wall.
- Weight capacity is 650 lbs. (295 kg).
- Call for freight quote.

| | Width x Length |
|-----------|----------------------|
| NC9203-47 | 4' x 7' (1.2 x 2.1m) |
| NC9203-57 | 5' x 7' (1.5 x 2.1m) |

CLINICA™ POWER PLATFORM

Foot pedal control for easy, hands-free use.

- Welded steel frame.
- Enclosed, single motor drive ensures quiet operation. 110 volt, 60Hz, 2.7 amps. Overload protection.
- Lifts from low to high in 15 seconds. Height adjusts from 19" to 28" (48 to 71cm).
- Premium, stain-resistant vinyl upholstery with knit-backing. Choice of eight colors (see page 20). Bumper guards protect vinyl.
- Lifting weight capacity is 600 lbs. (272 kg). Static weight capacity is 700 lbs. (317 kg).
- Call for freight quote.



| | Width x Length |
|-----------|----------------------|
| NC9211-47 | 4' x 7' (1.2 x 2.1m) |
| NC9211-57 | 5' x 7' (1.5 x 2.1m) |

Clinica™ Bariatrics Series Treatment Tables



Bariatrics Series tables are available in a variety of finishes and upholstery color options. See color and finish charts on page 20.

Both the hardwood six-leg and the all-welded steel base designs have a 600 lb. (272kg) weight capacity, these tables are built solid for heavy-duty use.

Bariatric Series Features

- Plywood tops attached with bolt through the top tee-nut fastening system.
 - 2" (5 cm) hi-density urethane foam padding.
 - Easy-to-clean laminate surfaces.
 - Solid hardwood legs.
 - Steel reinforced stress points.
 - Heavy-duty construction throughout.
 - 600 lb. (272kg) capacity.
 - Minimal assembly required.
- Power Table:**
- All-welded steel base and frame.



NEW! CLINICA™ BARIATRIC FULL CABINET TREATMENT TABLE

Three enclosed cabinets with adjustable shelves.

- Cabinet constructed with steel fasteners.
- Furniture-grade, laminated plywood cabinet frame. This model also available in gray.
- Six doors with concealed, self-closing Euro-style hinges.
- Measures 78" long x 31" high x 33" wide (198 x 79 x 84 cm).
- Three-year limited manufacturer warranty.
- Call for freight quote.

NC91019



NEW! CLINICA™ BARIATRIC POWER TABLE WITH POWER BACKREST

Heavy duty, power treatment table.

- All-welded steel base and frame.
- Dependable, single motor lift mechanism and motorized backrest with hands-free foot control.
- Easy-to-clean ABS plastic motor shroud.
- Hospital-grade plug. Overload protection. 110 volt/60 hz./3.7 amps.
- Measures 72" long x 20"-30" high x 33" wide (183 x 51-76 x 84 cm).
- One-year limited manufacturer warranty on frame and motor.
- Call for freight quote.

NC91017



NEW! CLINICA™ BARIATRIC AND LOW HEIGHT TREATMENT TABLES

Tables available in two different heights.

- Six-leg design with triple bolted corner legs.
- Built-in middle leg levelers.
- Steel reinforced frame.
- Laminate full length shelf with H-braces.
- Pillow included.
- Both tables measure 78" long x 33" wide (198 x 84 cm), heights are specific to each table.
- Three-year limited manufacturer warranty.
- Call for freight quote.

| | Height |
|-----------|-------------|
| NC91018-1 | 31" (79 cm) |
| NC91018-2 | 26" (66 cm) |

L-Series Features

- Solid hardwood legs (not laminated particle board).
- Option of Natural or Walnut Finish.
- Legs secured by H-braces and screws.
- Steel angle secures top to frame.
- Supports up to 350 lbs. (159kg).

S-Series Features

- Durable powder-coated steel frame.
- Supports up to 500 lbs. (227kg).
- Ships by UPS in two cartons. Minimal assembly required.
- Assembled units available. Call for details.

Clinica™ L-Series and S-Series Treatment Tables

L-Series Tables feature solid hardwood legs and steel connecting supports with a 350 lb. (159kg) weight capacity.

S-Series Tables have a powder-coated steel frame with a 500 lb. (227kg) weight capacity.

Both series feature an exclusive, flame-retardant foam that makes Clinica™ Treatment Tables stand apart from ordinary treatment tables. Available in a variety of table configurations, two table finishes and eight colors of germ, stain and scratch resistant premium vinyl upholstery (see page 20). All L-Series Treatment Tables have a three-year limited warranty against defects in workmanship. **Quick Ship** guarantees that tables will be shipped within three to five business days.

Premium upholstery



Solid hardwood legs



CLINICA™ MODEL 100

Standard treatment table with adjustable backrest.

- Included backrest adjusts from 0° to 74° in six fixed positions.
- Measures 72" long x 31" high (183 x 79 cm).
- S-Series tables are constructed of powder-coated steel and support up to 500 lbs. (227kg).
- L-Series tables support up to 350 lbs. (159kg).
- Available for Quick Ship.
- Call for freight quote.

| L-100 (Wood) | Width |
|---------------|-------------|
| NC9102-24 | 24" (61 cm) |
| NC9102-27 | 27" (69 cm) |
| NC9102-30 | 30" (76 cm) |
| S-100 (Steel) | Width |
| NC9107-27 | 27" (69 cm) |
| NC9107-30 | 30" (76 cm) |



Slate Blue upholstery and Natural finish shown.



Model L-100



Model S-100

S-Series basic steel treatment table.

CLINICA™ MODEL 200

Under-table shelf offers convenient storage space.

- Included backrest adjusts from 0° to 74° in six fixed positions.
- Measures 72" long x 31" high (183 x 79 cm).
- S-Series tables are constructed of powder-coated steel and support up to 500 lbs. (227kg).
- L-Series tables support up to 350 lbs. (159kg).
- Available for Quick Ship.
- Call for freight quote.

| L-200 (Wood) | Width |
|---------------|-------------|
| NC9103-24 | 24" (61 cm) |
| NC9103-27 | 27" (69 cm) |
| NC9103-30 | 30" (76 cm) |
| S-200 (Steel) | Width |
| NC9108-27 | 27" (69 cm) |
| NC9108-30 | 30" (76 cm) |



Slate Blue upholstery and Natural finish shown.



Model L-200



Model S-200

S-Series steel treatment table with shelf.

Finish and Upholstery Options

Clinica™ upholstered furniture is available in a variety of color selections. All Clinica™ heavy duty upholstery is coated with Permaplock3® for years of microbial, scratch and stain resistance. Choose from 8 attractive colors, shown below. Please indicate color ID number when placing order.


Upholstery color is required for all orders.

Table finish option is required for L-Series.

Clinica™ Upholstery Colors



Standard Base Laminate Colors for Frame, Legs, Shelving and Door



Gray is available as a third option for the Bariatric Full Cabinet Treatment Table (NC91019, pg. 18).

Concord Grape upholstery and Natural finish shown.



Model L-300

Model S-300

Model S-300D

S-Series steel treatment tables with tier shelves or drawers.

CLINICA™ MODEL 300

Open shelf and compartments provide easy access to storage.

- Included backrest adjusts from 0° to 74° in six fixed positions.
- Measures 72" long x 31" high (183 x 79 cm).
- S-Series tables are constructed of powder-coated steel and support up to 500 lbs. (227 kg)
- Available for Quick Ship.
- Call for freight quote.

| L-300 (Wood) | Width |
|----------------|-------------|
| NC9104-24 | 24" (61 cm) |
| NC9104-27 | 27" (69 cm) |
| NC9104-30 | 30" (76 cm) |
| S-300 (Steel) | Width |
| NC9109-27 | 27" (69 cm) |
| NC9109-30 | 30" (76 cm) |
| S-300D (Steel) | Width |
| NC9110-27 | 27" (69 cm) |
| NC9110-30 | 30" (76 cm) |

Concord Grape upholstery and Natural finish shown.



Model L-400

CLINICA™ MODEL 400

Convenient and discreet under-table storage.

- Included backrest adjusts from 0° to 74° in six fixed positions.
- Measures 72" long x 31" high (183 x 79 cm).
- Not available for Quick Ship.
- Call for freight quote.

| L-400 (Wood) | Width |
|--------------|-------------|
| NC9105-27 | 27" (69 cm) |
| NC9105-30 | 30" (76 cm) |



Model 100
Treatment Table

MEDALLION™ MODEL 100 TREATMENT TABLE

- Non-skid adjustable glides.
- One-piece top section measures 76" x 27" (193 x 69 cm).
- Height adjusts from 17" to 36" (43 to 91 cm).
- Call for freight quote.

NC9149



Model 200B
Treatment Table

MEDALLION™ MODEL 200B TREATMENT TABLE

- Two-piece top section measures 76" long x 34" wide (193 x 86 cm).
- Head section with contoured face opening measures 27" x 34" (69 x 86 cm). Raises 74°.
- One-piece body section measures 49" x 34" (134 x 86 cm).
- Height adjusts from 18" to 37" (46 to 94 cm).
- Maximum weight capacity is 500 lbs. (227 kg).
- Call for freight quote.

NC9154



Model 500R
Treatment Table

MEDALLION™ MODELS 500 AND 500R TREATMENT TABLES

- Five-piece top section measures 76" long x 27" wide (193 x 69 cm).
- Three-piece head section with contoured face opening measures 15" x 13" (38 x 33 cm). Raises 45° and lowers 90°.
- Armrests measure 12½" x 5½" (32 x 14 cm).
- Center section measures 22" x 27" (56 x 69 cm). Raises 30° with a spring loaded lock.
- Model 500 is the same as Model 500R, except the center section cannot be raised.
- Foot section measures 39" x 27" (99 x 69 cm).
- Call for freight quote.

NC9152 Model 500

NC9153 Model 500R

Medallion™ Treatment Tables—Strong, Stylish and Affordable.

Medallion™ Treatment Tables offer full height adjustment and section mobility so the patient can be placed in the proper position for maximum results with minimum effort. Adjustability allows greater patient comfort and ease getting on and off the table. Medallion™ Treatment Tables share these features:

- Heavy-duty vinyl coated Permacblock3® for years of microbial, scratch and stain resistance.
- Firm-density, foam tops measure 1½" (3.8 cm) thick.
- Radiused, seamless corners on upholstered components.
- Top sections have welded steel support frames.
- Pivoting joints and hardened steel shoulder bolts with locking nuts.
- Footswitch control.
- Power-assisted caster system (unless noted otherwise).
- 110 volts, 60Hz motor. (220 volts, 50Hz available on special order).
- Weight capacity is 400 lbs. (181 kg) (unless noted otherwise).

Medallion™ Treatment Tables Upholstery Color Options

Please specify color code when ordering treatment tables.

| | |
|------------------------|----------------------|
| Black 499-01 | Slate Blue 499-02 |
| Burgundy 499-03 | Gray 499-04 |
| Forest Green 499-05 | Navy 499-06 |
| Rose 499-07 | Taupe 499-08 |



Model 300R
Treatment Table

MEDALLION™ MODELS 300 AND 300R TREATMENT TABLES

- Three-piece top section measures 76" long x 27" wide (193 x 69 cm).
- One-piece head section with contoured face opening measures 15" x 27" (38 x 69 cm). Raises 45° and lowers 90°.
- Center section measures 22" x 27" (56 x 69 cm). Raises 30° with a spring-loaded lock.
- Model 300 is the same as Model 300R, except the center section cannot be raised.
- Foot section measures 39" x 27" (99 x 69 cm) and has 30° of movement.
- Call for freight quote.

NC9150 Model 300

NC9151 Model 300R



CLINICA™ ROUND HYDRAULIC LIFT TABLE

Tabletop can be raised, lowered, rotated and locked in place.

- Laminate top with tee-molded edge for added durability and protection.
- Top rotates and can be locked in place
- Tabletop measures 48" (122cm) in diameter.
- Hydraulic foot pump mechanism easily raises and lowers table top.
- Height adjusts from 27½" to 34" (70 to 86cm).
- Call for freight quote.

NC9133



CLINICA™ U-SHAPE TABLE

Treat up to four patients at a time.

- Laminate top with tee-molded edge for added durability and protection.
- Tubular steel legs adjust in height from 26" to 34" (66 to 86cm) in 1½" (3.8cm) increments. Positive-lock, spring plungers allow for easy height adjustments.
- Measures 18" (46cm) deep, 72" x 48" overall (183 x 122cm).
- Call for freight quote.

NC9132



CLINICA™ POWDER BOARD TABLES

Ideal for hand, arm and leg range-of-motion exercises.

- Smooth, scratch-resistant surface.
- Small table measures 40" x 29" x 7" (102 x 74 x 18cm).
- Large table measures 38" long x 63" (97 x 160cm) wide. Height adjusts from 26" to 34" (66 x 86cm) in 1½" (3.8cm) increments. Positive-lock, spring plungers allow for quick, easy height adjustments.
- Call for freight quote.



Small Table

NC9131 Large
NC9131-1 Small



CLINICA™ PORTABLE, HEIGHT ADJUSTABLE TABLE

Convenient, rolling, height-adjustable table accommodates most work areas.

- Gray laminate top with tee-molded edges for added protection and durability.
- Top measures 24" x 32" (61 x 81cm).
- Spring-assisted, dual columns adjust in height from 28" to 43" (71 to 109cm) in 1½" (3.8cm) increments.
- Chrome base features four casters.
- Weight capacity is 100 lbs. (45kg).
- Call for freight quote.

NC9130



CLINICA™ FIVE-LEG PNEUMATIC STOOL

Easy height adjustment with a single lever.

- Single-lever, pneumatic height adjustment from 19½" to 25" (50 to 64cm).
- Dual wheel, 2" (5.1cm), nylon casters.
- Seat measures 16" (41cm) in diameter and features 4" (10cm) thick padding.
- Available in eight easy-to-clean, stain-resistant vinyl colors (see page 20).
- BIFMA-tested, one-piece base is available in Black (377-21), Putty (377-22) or Slate (377-23).
- Call for freight quote.

NC9615

Foot Stool without Handle



FOOT STOOL

Steel-welded, chrome-plated construction.

Stool features non-skid rubber platform and anti-slip rubber tips for added security. Available with a handle for improved stability. Stool platform measures 14" wide x 1¼" deep x 9" high (36 x 3.2 x 23cm). Stool handle measures 36" (91cm) from floor.

NC93901-1 Without Handle

NC93901-2 With Handle

Foot Stool with Handle



CLINICA™ CHROME BASE STOOL

Height-adjustable stool with foot ring rolls easily from place to place.

- Easily spin the solid steel shaft to adjust height from 19" to 25" (48 to 64cm).
- Welded, 1" (2.5cm) tubular frame and foot ring are completely chrome-plated.
- Smooth-rolling casters measure 2" (5.1cm).
- Heavy-duty steel plate secures seat to base.
- Seat measures 14" (36cm) in diameter and features 4" (10cm) thick padding.
- Available in eight easy-to-clean, stain-resistant vinyl colors (see page 20).
- Call for freight quote.

NC9614



CLINICA™ PORTABLE MIRROR

Mobile mirror is ideal where space is limited.

- ANSI safety-backed mirror glass. Gray laminate construction.
- Measures 27" wide x 72" high (69 x 183cm). Casters measure 3" (7.6cm).
- Call for freight quote.

NC9500



CLINICA™ MOBILE MIRROR

Rolling mirror can be moved easily.


- ANSI safety-backed mirror glass. Natural, oak frame.
- Measures 27" wide x 72" high (69 x 183cm). Casters measure 3" (7.6cm).
- Call for freight quote.

NC9502

Clinica™ Parallel Bars Adjust Easily in Height and Width.




Telescoping steel tubing.



Spring-loaded plungers adjust height easily.

Two knobs at each upright ensure stability.



Adjust width easily on most models.

Clinica™ Parallel Bars feature quality construction with rounded, one-piece stainless steel handrails and bumper ends. Uprights are constructed of ultra-heavy, telescoping steel tube. Spring-loaded plungers with positive-lock positioning knobs provide easy height adjustment. A second knob at each upright ensures stability. Clinica™ Parallel Bars have a 350 lb. (159kg) weight capacity.

Some assembly required.



Folds flat for storage.

CLINICA™ WALL FOLDING PARALLEL BARS

An ideal choice where space is limited.

Handrails measure 7' (2.1m) long and fold up against the wall when not in use. Folded dimensions measure 14" high x 10" deep (36 x 25cm). Spring-loaded plungers adjust rail height easily from 28" to 42" (71 to 107cm) in 1½" (3.8cm) increments. Width at handrails is a fixed 25" (13cm). Call for freight quote.

| | Length |
|-----------|-----------|
| NC9304-07 | 7' (2.1m) |



Folding Parallel Bars

CLINICA™ FOLDING PARALLEL BARS

Easily adjustable system folds flat for storage.

Unique stabilizer bars stiffen the rails. Height adjusts from 22" to 35" (56 to 89cm) in 1½" (3.8cm) increments. Width adjusts from 16" to 24" (41 to 61cm) in 1½" (3.8cm) increments. Call for freight quote.

| | Length |
|-----------|------------|
| NC9303-07 | 7' (2.1m) |
| NC9303-10 | 10' (3.0m) |


NEW! BAILEY PLATFORM MOUNTED PARALLEL BARS

Handrails are mounted to a hardwood platform with easy-access ramps.

Solid birch plywood platform with durable laquer finish features 12" (30cm) long ramps for easy access. Chrome-plated steel bars are height adjustable from 26" to 44" (66 to 112cm) in 1" (2.5cm) increments. Fixed handrail width of 25½" (65cm). Available in a 10' or 12' (3.0 or 3.7m) handrail length. Platform measures 10' x 30" (3.0 x .76m). Weight capacity is 300 lbs. (136kg). The **Fixed Child's Handrail Attachment** is sold separately in 10' or 12' (3.0 or 3.7m) lengths and has 1" (2.5cm) diameter handrails for small hands. Call for freight quote on bars and handrail.

| Parallel Bars | Length |
|---------------|------------|
| NC93010-10 | 10' (3.0m) |
| NC93010-12 | 12' (3.7m) |

| Child's Handrail Attachment | Length |
|-----------------------------|------------|
| NC93011-10 | 10' (3.0m) |
| NC93011-12 | 12' (3.7m) |



Platform Mounted Parallel Bars

Fixed Child's Handrail Attachment



Floor Mount Parallel Bars

Width adjusts by 6" (15cm) on each side.

CLINICA™ FLOOR MOUNT PARALLEL BARS

Handrails install as a permanent fixture.

Height adjusts from 26" to 39" (66 to 99cm) in 1½" (3.8cm) increments. Call for freight quote.

| | Length |
|-----------|------------|
| NC9302-10 | 10' (3.0m) |
| NC9302-12 | 12' (3.7m) |

CLINICA™ TRAINING STAIRS WITH STEEL UPRIGHTS AND ROUNDED WOOD HANDRAILS

Ideal for ambulation and endurance training.

Available in three staircase styles: Small, One-Side and Convertible. All styles feature 9" (23 cm) deep steps. The **Small** style is a straight staircase with a 24" x 30" (61 x 76 cm) platform. One side features two steps with a 6" (15 cm) rise, while the other has three steps with a 4" (10 cm) rise. Footprint measures 36" x 51" (91 x 130 cm). The **One-Side** staircase has a 24" x 30" (61 x 76 cm) platform with three enclosed sides and four steps with a 6" (15 cm) rise. Footprint measures 36" x 60" (91 x 152 cm). The **Convertible** adjusts from a straight to a corner setup. Features a 30" x 36" (76 x 91 cm) platform. One side has four steps with a 6" (15 cm) rise while the other side has six steps with a 4" (10 cm) rise. Straight setup has a 45" x 104" (114 x 264 cm) footprint and the corner setup has a 47" x 58" (119 x 147 cm) footprint. All styles have a weight capacity of 350 lbs. (159 kg). Some assembly required. Call for freight quote.

NC9401 Small
NC9402 One-Side
NC9405 Convertible



Small



One-Side



Convertible:
Straight setup



Convertible:
Corner setup



CLINICA™ NESTING STOOLS

Multiple size stools nest for convenient storage.

These wooden stools are great for use at a desk, during seated mat activities or for improving body mechanics and posture while standing. Set of four stools includes: 10" x 12" x 2" (25 x 30 x 5.1 cm), 12" x 14" x 4" (30 x 36 x 10 cm), 14" x 18" x 6" (36 x 46 x 15 cm) and 16" x 20" x 8" (41 x 51 x 20 cm). Call for freight quote.

NC9431 Set of four

